## Alkaline phosphatase optimized

● Indicates analyzer(s) on which kit(s) can be used

Cat. no.	Bottle	Contents	704	717	736	747	902	904	911	914	917	MODULAR		
					737							P	D	
816388	1	Buffer, 6 x 42 ml												
	2	Buffer, 6 x 9 ml	●				●	●	●					
	2a	Substrate, 1 x 24 reagent-tablets												

Some analyzers and kits shown may not be available in all countries. For additional system applications, contact your local Roche Diagnostics representative.

**Intended use**

In vitro test for the quantitative determination of alkaline phosphatase (ALP) in human serum and plasma on automated clinical chemistry analyzers.

**Summary<sup>1-4</sup>**

Alkaline phosphatase in serum consists of four structural genotypes: the liver-bone-kidney type, the intestinal type, the placental type and the variant from germ cells. It occurs in osteoblasts, hepatocytes, the kidneys, spleen, placenta, prostate, leukocytes and the small intestine. The liver-bone-kidney type is particularly important.

A rise in the alkaline phosphatase activity occurs with all forms of cholestasis, particularly with obstructive jaundice. It is also elevated in diseases of the skeletal system, such as Paget's disease, hyperparathyroidism, rickets and osteomalacia, as well as with fractures and malignant tumors. A considerable rise in the alkaline phosphatase activity is sometimes seen in children and juveniles. It is caused by increased osteoblast activity following accelerated bone growth. Various reference values for the purposes of clinical evaluation have been assigned to differing age groups.

In 1946, Bessey, Lowry and Brock published a method for the determination of alkaline phosphatase using p-nitrophenyl phosphate as substrate buffered with glycine/NaOH. In 1967, Hausamen et al improved upon the method by using diethanolamine as buffer. The "optimized standard method" assay described here meets the 1972 recommendations of the "Deutschen Gesellschaft für Klinische Chemie (German Society of Clinical Chemistry)".

**Test principle<sup>4</sup>**

Colorimetric assay in accordance with a standardized method

- Sample and addition of R1 (buffer)
- Addition of R2 (buffer/substrate) and start of reaction:

In the presence of magnesium and zinc ions, p-nitrophenyl phosphate is hydrolyzed by phosphatases to form phosphate and p-nitrophenol. The p-nitrophenol released is proportional to the ALP activity and can be measured photometrically.

**Working solution concentration****R1** Buffer (bottle 1)

Diethanolamine buffer: 1.02 mol/l, pH 9.8; magnesium chloride: 0.51 mmol/l

**R2** Buffer/substrate (bottles 2 and 2a)

Diethanolamine buffer: 1.02 mol/l, pH 9.8; magnesium chloride: 0.51 mmol/l; p-nitrophenyl phosphate: 61 mmol/l

**Precautions and warnings**

For in vitro diagnostic use.

Exercise the normal precautions required for handling all laboratory reagents.

**Reagent handling**

R1: Ready for use.

R2: To prepare 9 ml, add 4 reagent tablets from bottle 2a to one bottle 2. Dissolve by swirling gently.

**Storage and stability**

Unopened kit components: Up to the expiration date at 2–8°C

R1: 14 days opened and refrigerated on the analyzer

R2: 14 days opened and refrigerated on the analyzer.

**Specimen collection and preparation**

Collect serum using standard sampling tubes.

Heparinized plasma

Stability<sup>5</sup>: 2 days at 20–25°C

7 days at 4–8°C

4 weeks at –20°C

Centrifuge samples containing precipitate before performing the assay.

**Testing procedure***Materials provided*

- Working solutions as described above

*Additional materials required*

- Calibrators and controls as indicated below
- 0.9% NaCl
- Chimneys

**Assay**

Refer to the appropriate operator's manual and/or Instrument Settings section of this insert for analyzer-specific assay instructions. The performance of applications not validated by Roche is not warranted and must be defined by the user.

**Calibration**

Absorption of atmospheric CO<sub>2</sub> by the opened reagent bottle leads to impaired calibration stability. This kit therefore requires the use of color-coded chimneys, which reduce the uptake of CO<sub>2</sub> by the reagents. The chimneys should be placed directly into the appropriate reagents: white for R1, black for R2. The chimneys can be reused for reagent bottles within the same kit.

S1: 0.9% NaCl

S2: C.f.a.s. (Calibrator for automated systems)

Determine the K-factor when setting up the instrument.

Calibration frequency

Two-point calibration is recommended:

- as required following quality control procedures

Calibration verification: Not necessary.

**Quality control**

For quality control, use Precinorm U, Precinorm E, Precipath U, Precipath E or other suitable control material. The control intervals and limits should be adapted to the individual laboratory and country-specific requirements. Values obtained should fall within the established limits. Each laboratory should establish corrective measures to be taken if values fall outside the limits.

**Calculation**

The Roche/Hitachi systems automatically calculate the alkaline phosphatase activity of each sample.

Conversion factor: U/l x 0.0167 = μkat/l

**Limitation – interference<sup>6,7</sup>**

Criterion: Recovery within ± 10% of initial value.

Icterus: No significant interference up to an I index of 60 (approximate conjugated and unconjugated bilirubin concentration: 60 mg/dl).

Hemolysis: No significant interference up to an H index of 1000 (approximate hemoglobin concentration: 1000 mg/dl).

Lipemia (Intralipid): No significant interference up to an L index of 1000 (approximate triglycerides concentration: 2000 mg/dl). There is poor correlation between turbidity and triglycerides concentration.

Close bottle 2a well after removing the tablets.

Close bottle 2a well after removing the tablets.

Close bottle 2a well after removing the tablets.

Close bottle 2a well after removing the tablets.

**Measuring range**

Measuring range: 5–2000 U/l (0.08–33.34 μkat/l)

Determine samples having higher activities via the rerun function. On instruments without rerun function, manually dilute the samples with 0.9% NaCl or distilled/deionized water (e.g. 1 + 4). Multiply the result by the appropriate dilution factor (e.g. factor 5).

## Reference range according to Rosalki

		25°C <sup>§</sup>	30°C*	37°C*
U/l	Men	< 180	< 220	< 270
	Women	< 160	< 195	< 240
µkat/l	Men	< 3.00	< 3.67	< 4.50
	Women	< 2.67	< 3.25	< 4.00

\*Calculated values

The following factors were used for the conversion of reference values determined at 25°C: 1.24 (30°C) and 1.52 (37°C).

For reference values for children, please refer to "Reference Ranges for Adults & Children, Preanalytical Values". Heil W, Koberstein R, Zawta B (published by Boehringer Mannheim GmbH, 1997).

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference range. For diagnostic purposes, the alkaline phosphatase results should always be assessed in conjunction with the patient's medical history, clinical examination and other findings.

### Specific performance data

The data determined using a Roche/Hitachi system are given below. Results obtained in individual laboratories may differ.

### Imprecision<sup>§</sup>

Reproducibility was determined using human samples and controls in an internal protocol (within run n = 21, between day n = 10). The following results were obtained.

Sample	Within run			Between day		
	Mean U/l	SD U/l	%CV	Mean U/l	SD U/l	%CV
Human serum	458	2.23	0.5	357	7.80	2.2
Precinorm U	457	1.66	0.4	439	10.10	2.3
Precipath U	579	2.39	0.4	563	12.10	2.1

### Analytical sensitivity (lower detection limit)<sup>§</sup>

5 U/l or 0.08 µkat/l

The detection limit represents the lowest measurable alkaline phosphatase activity that can be distinguished from zero. It is calculated as three standard deviations of 21 replicates of the lowest standard.

### Method comparison<sup>§</sup>

A comparison of the new wavelength combination 546/450 nm (y) and the old combination 700/415 nm (x) using the Roche alkaline phosphatase opt. reagents on the Roche/Hitachi 911 gave the following correlation (U/l):

Passing/Bablok <sup>9,10</sup>	Linear regression
$y = 7.369 + 1.012 x$	$y = 17.445 + 0.994 x$
$r = 0.9997$	$r = 0.9997$
SD (md 95) = 15.2	$Sy.x = 6.6$
Number of samples measured: 60	

The activities of the samples were between 213 and 1820 U/l.

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## Instrument settings

### Roche/Hitachi 904/911/912

Read in the application parameters from the application diskette or barcode sheet, as appropriate.

### Roche/Hitachi 704

Temperature: 25°/30°/37°C

PROGRAM 2 CHEMISTRY PARAMETERS	
TEST	[ALP]
ASSAY CODE	[2(RATE-A)]-[19]-[32]
SAMPLE VOLUME	[9]
R1 VOLUME	[350]-[50]-[NO]
R2 VOLUME	[70]-[20]-[NO]
WAVELENGTH	[546]-[450]
CALIB. TYPE	[LINEAR]-[0]-[0]
STD. (1) CONC.-POS.	[ ]-[ ]
STD. (2) CONC.-POS.	[ ]-[ ]
STD. (3) CONC.-POS.	[0]-[0]
STD. (4) CONC.-POS.	[0]-[0]
STD. (5) CONC.-POS.	[0]-[0]
STD. (6) CONC.-POS.	[0]-[0]
UNIT	[ ]
SD LIMIT	[0.1]
DUPLICATE LIMIT	[100]
SENSITIVITY LIMIT	[0]
ABS. LIMIT (INC/DEC)	[14000]-[INCREASE]
PROZONE LIMIT	[0]-[LOWER]
NORMAL RANGE	[ ]-[ ]
INSTRUMENT FACTOR	[1.00]

— Data entered by the operator

## Roche/Hitachi 902

No.	<Chemistry>	
1	Test Name	ALP
2	Assay Code(Mthd)	RATE A
3	Assay Code (2. Test)	0
4	Reaction Time	10
5	Assay Point 1	26
6	Assay Point 2	35
7	Assay Point 3	0
8	Assay Point 4	0
9	Wavelength (SUB)	546
10	Wavelength (MAIN)	450
11	Sample Volume	7.0
12	R1 Volume	250
13	R1 Pos.	.....
14	R1 Bottle Size	Large
15	R2 Volume	0
16	R2 Pos.	0
17	R2 Bottle Size	Small
18	R3 Volume	50
19	R3 Pos.	.....
20	R3 Bottle Size	Small
21	Calibration Type (Type)	Linear
22	Calibration Type (Wght)	0
23	Calib. Conc. 1	0
24	Calib. Pos. 1	.....
25	Calib. Conc. 2	.....
26	Calib. Pos. 2	.....
27	Calib. Conc. 3	0
28	Calib. Pos. 3	0
29	Calib. Conc. 4	0
30	Calib. Pos. 4	0
31	Calib. Conc. 5	0
32	Calib. Pos. 5	0
33	Calib. Conc. 6	0
34	Calib. Pos. 6	0
35	S1 ABS	0
36	K Factor	10000
37	K2 Factor	10000
38	K3 Factor	10000
39	K4 Factor	10000
40	K5 Factor	10000
41	A Factor	0
42	B Factor	0
43	C Factor	0
44	SD Limit	0.1
45	Duplicate Limit	30
46	Sens. Limit	240
47	S1 ABS. Limit (L)	-32000
48	S1 ABS. Limit (H)	32000
49	ABS. Limit	14000
50	ABS. Limit (D/I)	Increase
51	Prozone Limit	0
52	Proz Limit (Upp/Low)	Lower
53	Prozone (End Point)	35
54	Expect. Value (L)	.....
55	Expect. Value (H)	.....
56	Instr. Fact. (a)	1.0
57	Instr. Fact. (b)	0.0
58	Key Setting	.....

.... Data entered by the operator

For detailed information, consult the operator manuals for Roche/Hitachi systems, the respective application sheets and the package inserts for the calibrators and control sera.

Precinorm and Precipath are trademarks of a member of the Roche Group.

Intralipid is a trademark of KabiPharmacia, Inc.

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© 1999 Roche Diagnostics
Roche Diagnostics GmbH, D-68298 Mannheim, Germania
Roche Diagnostics Corporation, Indianapolis, IN, USA