## Bilirubin DPD

● Indicates Roche/Hitachi analyzer(s) on which kit(s) can be used

Cat. No.	Bottle	Contents	704	717	737	747	902	904	911 912	914	917	MODULAR P	D
1552414	1	Detergent/hydrochloric acid, 9 x 62 ml											
	2	Diazo reagent, 9 bottles of granulate for 15.5 ml each									●		
1489194	1	Detergent/hydrochloric acid, 12 x 50 ml	●				●	●	●				
	2	Diazo reagent, 6 bottles of granulate for 22 ml each											
1489429	1	Detergent/hydrochloric acid, 6 x 100 ml		●				●	●				
	2	Diazo reagent, 3 bottles of granulate for 46 ml each											

Some analyzers and kits shown may not be available in all countries. For additional system applications, contact your local Roche representative.

### Intended use

For the in vitro quantitative determination of total bilirubin in serum and plasma of adults and newborn on automated clinical chemistry analyzers.

### Summary<sup>1-4</sup>

Bilirubin is produced during normal and abnormal degradation of erythrocytes in the reticuloendothelial system. Bilirubin determinations are made use of in the diagnosis of liver diseases, the detection of hemolytic anemia and in the assessment of the degree of severity of jaundice.

Since the introduction of the diazo method for determining bilirubin by Ehrlich in 1883<sup>1</sup>, a number of modifications for improving the reaction have been proposed. In the Malloy-Evelyn method<sup>2</sup>, methanol is used as a catalyst for the azo-coupling of indirect bilirubin. In addition it keeps azobilirubin in solution. A serious disadvantage of this method is that protein can be precipitated by the methanol, thus possibly leading to falsely depressed values.

In 1938, Jendrassik and Grof<sup>3</sup> presented a test yielding reliable results. However, this method is time-consuming and involves a number of pipetting steps.

The determination presented here corresponds to the method developed by Wahlefeld et al<sup>4</sup> in which a detergent is used to accelerate the reaction and prevent the precipitation of protein. The diazo reagent used is 2,5-dichlorophenyl diazonium tetrafluoroborate (DPD), which under acidic conditions couples very rapidly with bilirubin. The procedure is simple and correlates well with the Jendrassik-Grof method.

### Test principle<sup>4</sup>

Colorimetric assay

- Sample and addition of R1 (detergent/hydrochloric acid)
- Addition of R2 (diazo reagent) and start of the reaction:  
Indirect bilirubin is liberated by the detergent. In strong acid solution containing 2,5-dichlorophenyl diazonium salt, total bilirubin couples to form azobilirubin.

The color intensity of the red azo dye formed is directly proportional to the total bilirubin and can be determined photometrically.

### Working solution concentration

**R1** Detergent/hydrochloric acid

Detergent; hydrochloric acid: 120 mmol/l

**R2** Diazo reagent

2,5-dichlorophenyl diazonium salt: 3.0 mmol/l

### Precautions and warnings

For in vitro diagnostic use.

Exercise the normal precautions required for handling all laboratory reagents.

**CAUTION. WARNING.** Bottle 1 contains hydrochloric acid which can cause caustic burns. Avoid contact with eyes, skin and mucous membranes. Flush affected areas with copious amounts of water. Get immediate medical attention for eyes, or if inhaled or ingested.

### Reagent handling

1552414/1489194/1489429

R1: Contents are ready for use

1552414/1489194

R2: Add approximately 10 ml of dist. water to one bottle 2 and dissolve the granulate by swirling the bottle gently. Make up to the mark with distilled water and mix. Do not stopper the bottle containing diazo reagent solution with a desiccant cap.

1489429

R2: Add approximately 30 ml of dist. water to one bottle 2 and dissolve the granulate by swirling the bottle gently. Make up to the mark with distilled water and mix. Do not stopper the bottle containing diazo reagent solution with a desiccant cap.

### Storage and stability

Unopened kit components: Up to the expiration date at 2–8°C

R1: 42 days opened and refrigerated on the analyzer

R2: 14 days opened and refrigerated on the analyzer

### Specimen collection and preparation

Collect serum using standard sampling tubes.

Heparin or EDTA plasma

Stability<sup>5</sup>: Store samples protected from light. Perform the assay immediately.

Centrifuge samples containing precipitate before performing the assay.

### Testing procedure

*Materials supplied*

- Working solution as described above.

*Additional materials required*

- Calibrators and controls as indicated below
- 0.9% NaCl

### Assay

Refer to the appropriate operator's manual and/or the Instrument Settings section of this insert for analyzer specific assay instructions. The performance of applications not validated by Roche is not warranted and must be defined by the user.

### Calibration

Standardization<sup>6</sup>: The bilirubin method has been calibrated against the "Dourmas candidate reference method"/SRM 916a.

S1: 0.9% NaCl

S2: C.f.a.s. (Calibrator for automated systems)

Calibration frequency

Recalibration is recommended:

- as a blank calibration after 24 hours
- as a blank calibration after reagent bottle change
- as a two point calibration after reagent lot change
- as a two point calibration if required following quality control procedures

Calibration verification: Not necessary


### Quality control

For quality control, use Precinorm U, Precinorm U plus, Precipath U, Precipath U plus, Precibil or other suitable control material. The control intervals and limits should be adapted to the individual laboratory and country-specific requirements. Values obtained should fall within established limits. Each laboratory should establish corrective measures to be taken if values fall outside the limits.

### Calculation

Roche/Hitachi systems automatically calculate the total bilirubin concentration of each sample.

Conversion factor: mg/dl x 17.1 = μmol/l

## Limitations – interference<sup>10,11</sup>

Occasionally, somewhat higher values for direct bilirubin than for total bilirubin may be yielded by samples in which the total and direct bilirubin concentrations are similar. In such samples, the bilirubin concentration determined by both methods should be stated.

Criterion: Recovery within  $\pm 10\%$  of initial value.

Lipemia interferes.

Hemolysis interferes in samples from adults. Hemolysis (up to approx. 350 mg/dl hemoglobin) does not cause interference in neonatal samples.<sup>8,9</sup>

Indican: No significant interference by indican up to concentrations of 0.2 mmol/l. At higher indican concentrations, which can occur in renal failure, spuriously elevated bilirubin values are measured.

The effect of 18 commonly used pharmaceuticals was tested in vitro. Of these, only cyclosporine at therapeutic concentrations led to lower recovery figures (values approximately 10% depressed).

## Measuring/reportable range

0.1–30 mg/dl or 1.71–315  $\mu\text{mol/l}$

Determine samples with higher concentrations via the rerun function. On instruments without rerun function, manually dilute samples with higher concentrations using 0.9% NaCl or distilled/deionized water (e.g. 1 + 3). Multiply the result by the appropriate factor (e.g. 4).

## Reference range

Reference range according to Thomas

Total bilirubin<sup>10</sup>: up to 1.1 mg/dl or 18.8  $\mu\text{mol/l}$

Reference range according to Sherlock and Meites

Adults and children<sup>11</sup>: up to 1.0 mg/dl or 17.1  $\mu\text{mol/l}$

Newborn<sup>12</sup>

Age of newborn	Premature
24 hours	1.0–6.0 mg/dl or 17.1–102.6 $\mu\text{mol/l}$
48 hours	6.0–8.0 mg/dl or 102.6–136.8 $\mu\text{mol/l}$
3–5 days	10.0–15.0 mg/dl or 171–256.5 $\mu\text{mol/l}$
Age of newborn	Full term
24 hours	2.0–6.0 mg/dl or 34.2–102.6 $\mu\text{mol/l}$
48 hours	6.0–7.0 mg/dl or 102.6–119.7 $\mu\text{mol/l}$
3–5 Tage	4.0–12.0 mg/dl or 68.4–205.2 $\mu\text{mol/l}$

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference range. For diagnostic purposes, total bilirubin results should always be assessed in conjunction with the patient's medical history, clinical examinations and other findings.

## Specific performance data

The data determined using a Roche/Hitachi system are given below. Results obtained in individual laboratories may differ.

## Imprecision<sup>6</sup>

Reproducibility was determined using human samples and controls in an internal protocol (n = 21). The following results were obtained.

Sample	Within run			Between day		
	Mean mg/dl	SD mg/dl	%CV	Mean mg/dl	SD mg/dl	%CV
Human serum	2.1	0.026	1.3	2.1	0.041	1.9
Precinorm U	2.1	0.038	1.8	2.1	0.044	2.1
Precipath U	5.0	0.060	1.2	5.0	0.079	1.6

## Analytical sensitivity<sup>6</sup>

0.1 mg/dl or 1.71  $\mu\text{mol/l}$

The lower detection limit represents the lowest measurable total bilirubin concentration that can be distinguished from zero. It is calculated as three standard deviations of 21 replicates of the lowest standard.

## Method comparison<sup>6</sup>

A comparison of the bilirubin assay using Roche Bilirubin DPD granulate reagent (y) with Roche Bilirubin DPD lyophilized reagent (x) on a Roche/Hitachi 704 analyzer gave the following correlation (mg/dl):

Passing/Bablok <sup>13,14</sup>	Linear regression
$y = 0.017 + 1.023 x$	$y = -0.009 + 1.031 x$
$r = 1.000$	$r = 1.000$
SD (md 95) = 0.104	Sy.x = 0.054

Number of samples measured: 50


The activities of the samples were between 0.19 and 14.29 mg/dl.

## Please note

The determination can also be performed at 25°C or 30°C.

## References

- Ehrlich P. Charite Ann 1883;8:140.
- Malloy HT, Evelyn KA. The determination of bilirubin with the photoelectric colorimeter. J Biol Chem 1937;119:481–490.
- Jendrassik L, Grof P. Vereinfachte photometrische Methoden zur Bestimmung des Bilirubins. Biochem Z 1938;297:81–89.
- Wahlefeld AW, Herz G, Bernt E. Modification of the Malloy-Evelyn method for a simple, reliable determination of total bilirubin in serum. Scand J clin Lab Invest 1972;29 Supplement 126:Abstract 11.12.
- Tietz NW, ed. Clinical Guide to Laboratory Tests, 3rd ed. Philadelphia, PA: WB Saunders, 1995:88.
- Data on file.
- Glick MR, Ryder KW, Jackson SA. Graphical Comparisons of Interferences in Clinical Chemistry Instrumentation. Clin Chem 1986;32:470–474.
- Mazzachi BC, Peake MJ. How reliable is the Hitachi DPD bilirubin method? Poster Presentations. Clin Biochem Revs 1995;16:73.
- Schlebusch H, Schneider C, Liappis N et al. Bilirubin determinations in neonatal sera: precision, accuracy and sensitivity to hemoglobin interference of six routine methods. AACC Poster Abstract. Clin Chem 1995;41 Supplement S6:S95.
- Thomas H, ed. Labor und Diagnose, 4th ed. Marburg: Die Medizinische Verlagsgesellschaft 1992:235.
- Sherlock S. Liver Disease. London: Churchill, 1951.
- Meites S. Pediatric Clinical Chemistry: A Survey of Normals, Methods, and Instrumentation, with Commentary. Washington, DC: 1977:48.
- Passing H, Bablok W. A New Biometrical Procedure for Testing the Equality of Measurements from Two Different Analytical Methods. J Clin Chem Clin Biochem 1983;21:709–720.
- Bablok W et al. A General Regression Procedure for Method Transformation. J Clin Chem Clin Biochem 1988;26:783–790.

## Instrument settings

### Roche/Hitachi 917

Read in the application parameters from the application diskette or barcode sheet, as appropriate.

### Roche/Hitachi 704

Temperature: 37°C

PROGRAM 2 CHEMISTRY PARAMETERS	
TEST	[BIL-T]
ASSAY CODE	[2(2POINT)]-[15]-[32]
SAMPLE VOLUME	[10]
R1 VOLUME	[350]-[50]-[NO]
R2 VOLUME	[70]-[20]-[NO]
WAVELENGTH	[660]-[570]
CALIB. METHOD	[LINEAR]-[0]-[0]
STD. (1) CONC.-POS.	[ ]-[ ]
STD. (2) CONC.-POS.	[ ]-[ ]
STD. (3) CONC.-POS.	[0]-[0]
STD. (4) CONC.-POS.	[0]-[0]
STD. (5) CONC.-POS.	[0]-[0]
STD. (6) CONC.-POS.	[0]-[0]
UNIT	[ ]
SD LIMIT	[0.1]
DUPLICATE LIMIT	[100]
SENSITIVITY LIMIT	[0]
ABS. LIMIT (INC/DEC)	[0]-[INCREASE]
PROZONE LIMIT	[32000]-[UPPER]
EXPECTED VALUE	[ ]-[ ]
INSTRUMENT FACTOR	[1.00]

\_\_\_ Data entered by the operator

### Roche/Hitachi 717

Temperature: 37°C

PROGRAM 2 CHEMISTRY PARAMETERS	
TEST	[BIL-T]
ASSAY CODE	[2(2POINT)]-[24]-[50]
SAMPLE VOLUME	[7]-[3]
R1 VOLUME	[250]-[100]-[NO]
R2 VOLUME	[50]-[50]-[NO]
WAVELENGTH	[660]-[570]
CALIB. METHOD	[LINEAR]-[0]-[0]
STD. (1) CONC.-POS.	[ ]-[ ]
STD. (2) CONC.-POS.	[ ]-[ ]
STD. (3) CONC.-POS.	[0]-[0]
STD. (4) CONC.-POS.	[0]-[0]
STD. (5) CONC.-POS.	[0]-[0]
STD. (6) CONC.-POS.	[0]-[0]
SD LIMIT	[0.1]
DUPLICATE LIMIT	[100]
SENSITIVITY LIMIT	[0]
ABS. LIMIT (INC/DEC)	[0]-[INCREASE]
PROZONE LIMIT	[32000]-[UPPER]
EXPECTED VALUE	[ ]-[ ]
PANIC VALUE	[ ]-[ ]
INSTRUMENT FACTOR	[1.00]

\_\_\_ Data entered by the operator

### Roche/Hitachi 902

No.	<Chemistry>	
1	Test Name	BIL-T
2	Assay Code (Mthd)	2 Point End
3	Assay Code (2. Test)	0
4	Reaction Time	10
5	Assay Point 1	17
6	Assay Point 2	35
7	Assay Point 3	0
8	Assay Point 4	0
9	Wavelength (SUB)	700
10	Wavelength (MAIN)	570
11	Sample Volume	7.0
12	R1 Volume	250
13	R1 Pos.	.....
14	R1 Bottle Size	Large
15	R2 Volume	0
16	R2 Pos.	0
17	R2 Bottle Size	Small
18	R3 Volume	50
19	R3 Pos.	.....
20	R3 Bottle Size	Small
21	Calib. Type (Type)	Linear
22	Calib. Type (Wght)	0
23	Calib. Conc. 1	0.00
24	Calib. Pos. 1	.....
25	Calib. Conc. 2	.....
26	Calib. Pos. 2	.....
27	Calib. Conc. 3	0
28	Calib. Pos. 3	0
29	Calib. Conc. 4	0
30	Calib. Pos. 4	0
31	Calib. Conc. 5	0
32	Calib. Pos. 5	0
33	Calib. Conc. 6	0
34	Calib. Pos. 6	0
35	S1 ABS.	0
36	K Factor	10000
37	K2 Factor	10000
38	K3 Factor	10000
39	K4 Factor	10000
40	K 5 Factor	10000
41	A Factor	0
42	B Factor	0
43	C Factor	0
44	SD Limit	0.1
45	Duplicate Limit	35
46	Sens. Limit	240
47	S1ABS. Limit (L)	-32000
48	S1ABS. Limit (H)	32000
49	ABS. Limit	0
50	ABS. Limit (D/I)	Increase
51	Prozone Limit	32000
52	Proz Limit (Upp/Low)	Upper
53	Prozone (End Point)	35
54	Expect. Value (L)	.....
55	Expect. Value (H)	.....
56	Instr. Fact. (a)	1.0
57	Instr. Fact. (b)	0.0
58	Key Setting	.....

.... Data entered by the operator

For detailed information, consult the operator manuals for Roche/Hitachi systems, the respective application sheets and the package inserts for the calibrators and control sera. Precinorm, Precipath and Precibil are trademarks of a member of the Roche Group.  
©1999 Roche Diagnostics

Roche Diagnostics GmbH, D-68298 Mannheim, Germany  
Roche Diagnostics Corporation, Indianapolis, IN, USA  
US Customer Technical Support 1-800-428-2336

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
## Bibliografia

- Ehrlich P. Charite Ann 1883;8:140.
- Malloy HT, Evelyn KA. The determination of bilirubin with the photoelectric colorimeter. J Biol Chem 1937;119:481–490.
- Jendrossik L, Grof P. Vereinfachte photometrische Methoden zur Bestimmung des Bilirubins. Biochem Z 1938;297:81–89.
- Wahlefeld AW, Herz G, Bernt E. Modification of the Malloy-Evelyn method for a simple, reliable determination of total bilirubin in serum. Scand J Clin Lab Invest 1972;29 Supplement 126:Ab-stract 11.12.
- Tietz NW (Ed.). Clinical Guide to Laboratory Tests, 3ª edição, Philadelphia, PA: WB Saunders, 1995:88.
- Documentação da Roche.
- Glick MR, Ryder KW, Jackson SA. Graphical Comparisons of Interferences in Clinical Chemistry Instrumentation. Clin Chem 1986;32:470–474.
- Mazzachi BC, Peake MJ. How reliable is the Hitachi DPD bilirubin method? Poster Presentations. Clin Biochem Revs 1995;16:73.
- Schlebusch H, Schneider C, Liappis N et al. Bilirubin determinations in neonatal sera: precision, accuracy and sensitivity to hemoglobin interference of six routine methods. AACC Poster Abstract. Clin Chem 1995;41 Supplement S6:S95.
- Thomas H (Ed.). Labor und Diagnose, 4. Auflage. Marburg: Die Medizinische Verlagsgesellschaft 1992:235.
- Sherlock S. Liver Disease. London: Churchill, 1951.
- Meites S. Pediatric Clinical Chemistry: A Survey of Normals, Methods, and Instrumentation, with Commentary. Washington, DC: 1977:48.
- Passing H, Bablok W. A New Biometrical Procedure for Testing the Equality of Measurements from Two Different Analytical Methods. J Clin Chem Clin Biochem 1983;21:709–720.
- Bablok W et al. A General Regression Procedure for Method Transformation. J Clin Chem Clin Biochem 1988;26:783–790.

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Roche Diagnostics GmbH, D-68298 Mannheim, Alemanha