Tina-quant® [a] CRP

● Indicates Roche/Hitachi analyzer(s) on which kit(s) can be used

Cat. no.	Bottle	Contents	704	717	747	747-400	902	904	911 912	914	917	MODULAR	
												P	D
1551922	1	Buffer, 6 x 20 ml										●	●
	2	Anti-CRP antibody/buffer, 6 x 4 ml											
1730371	1	Buffer, 6 x 66 ml										●	●
	2	Anti-CRP antibody/buffer, 6 x 11 ml											
1929372	1	Buffer, 4 x 267 ml											●
	2	Anti-CRP antibody/buffer, 4 x 57 ml											●
1299859	1	Buffer, 6 x 16 ml	●	●			●	●	●	●			
	2	Anti-CRP antibody/buffer, 1 x 16 ml											
1489054	1	Buffer, 6 x 100 ml		●				●	●	●			
	2	Anti-CRP antibody/buffer, 2 x 50 ml											
1776371	1	Buffer, 5 x 300 ml			●	●							
1776428	2	Anti-CRP antibody/buffer, 2 x 180 ml											

Some analyzers and kits shown may not be available in all countries. For additional system applications, contact your local Roche representative.

Intended use

Immunoturbidimetric assay for the in vitro quantitative determination of CRP in human serum and plasma on automated clinical chemistry analyzers.

Summary¹⁻⁵

C-reactive protein is the classic acute phase protein to inflammatory reactions. It is synthesized by the liver and consists of five identical polypeptide chains forming a five-membered ring of molar mass 120000 daltons. CRP is the most sensitive of the acute phase reactants, and its concentration increases rapidly during inflammatory processes. Complexed CRP activates the complement system beginning with C1q. CRP then initiates opsonization and phagocytosis of invading cells, but its main function is to bind and detoxify endogenous toxic substances produced as a result of tissue damage.


CRP assays are used to detect systemic inflammatory processes (apart from certain types of inflammation such as SLE and Colitis ulcerosa); to assess treatment of bacterial infections with antibiotics; to detect intrauterine infections with concomitant premature amniorrhexis; to differentiate between active and inactive forms of disease with concurrent infection, e.g. in patients suffering from SLE or Colitis ulcerosa; to therapeutically monitor rheumatic diseases and assess anti-inflammatory therapy; to determine the presence of postoperative complications at an early stage, such as infected wounds, thrombosis and pneumonia; and to distinguish between infection and bone-marrow transplant rejection.

Various assay methods are available for CRP determination such as nephelometry and turbidimetry. The Roche CRP assay is based on the principle of immunological agglutination.

Test principle⁶⁻⁸

Immunoturbidimetric assay

- Sample and addition of R1 (buffer)
- Addition of R2 (anti-CRP antibody/buffer) and start of reaction:

Anti-CRP antibodies react with antigen in the sample to form an antigen/antibody complex. Following agglutination, this is measured turbidimetrically. Addition of PEG allows the reaction to progress rapidly to the end point, increases sensitivity and reduces the risk of samples containing excess antigen producing false negative results.

Working solution concentration

R1 Buffer

TRIS/HCl buffer*: 100 mmol/l, pH 7.5; NaCl: 300 mmol/l; polyethylene glycol: 2.0%; preservative

R2 Anti-CRP antibody/buffer

Anti-human anti-CRP antibody (goat): dependent on titer; TRIS/HCl buffer*: 100 mmol/l, pH 8.0; NaCl: 300 mmol/l; polyethylene glycol: 2.0%; preservative

*TRIS = Tris(hydroxymethyl)-aminomethane hydrochloride

Precautions and warnings

For in vitro diagnostic use.

Exercise the normal precautions required for handling all laboratory reagents.

- Disposal of all waste material should be in accordance with local guidelines.

Reagent handling

R1: Ready for use

R2: Ready for use

Storage and stability

Unopened kit components: Up to the expiration date at 2–8°C

R1: 90 days opened and refrigerated on the analyzer

R2: 90 days opened and refrigerated on the analyzer

Specimen collection and preparation

Collect serum using standard sampling tubes.

Li-heparin or EDTA plasma

Stability⁹: 3 days at 20–25°C

8 days at 4–8°C

3 years at -20°C

Centrifuge samples containing precipitate before performing the assay.

Testing procedure

Materials provided

- Working solutions as described above

Additional materials required

- Calibrators and controls as indicated below
- 0.9% NaCl
- General laboratory equipment

Assay

Refer to the appropriate operator's manual and/or the Instrument Settings section of this insert for analyzer specific assay instructions. The performance of applications not validated by Roche is not warranted and must be defined by the user.

Calibration

Standardization: the CRP method was standardized against CRM 470 (RPPHS - Reference Preparation for Proteins in Human Serum).

Roche/Hitachi 704/717/747/902/914/MODULAR D

S1: 0.9% NaCl

S2–6: Preciset Serum proteins, Cat. No. 1876406

Roche/Hitachi 904/911/912/917/MODULAR P

S1: 0.9 % NaCl

S2: C.f.a.s. Proteins (Calibrator for automated systems), Cat. No. 1355279.

Multiply the lot-specific C.f.a.s. Proteins calibrator value by the factors below to determine the standard concentrations for the six-point calibration curve:

S2: 0.05 S5: 1.00

S3: 0.10 S6: 2.50

S4: 0.20

Calibration frequency

Full calibration is recommended:

- after reagent lot change
 - as required following quality control procedures
- Calibration verification: Not necessary

Quality control

For quality control use Precinorm Protein, Precipath Protein or other suitable control material. The control intervals and limits must be adapted to the individual laboratory requirements. Values obtained should fall within established limits. Each laboratory should establish corrective measures to be taken if values fall outside the limits.

Calculation

Roche/Hitachi systems automatically calculate the CRP concentration of each sample.

Conversion factor: mg/dl x 0.01 = g/l mg/dl x 10 = mg/l
g/l x 100 = mg/dl mg/l x 0.1 = mg/dl

Limitations – interference¹⁰

Criterion: Recovery within ± 10% of initial values.

Icterus: No significant interference up to an I index of 60 (approximate conjugated and unconjugated bilirubin concentration: 60 mg/dl or 1026 µmol/l).

Hemolysis: No significant interference up to an H index of 700. (approximate hemoglobin concentration: 700 mg/dl or 435 µmol/l)

Lipemia (Intralipid): No significant interference up to an L index of 1000 (approximate triglyceride concentration: 2000 mg/dl or 22.8 mmol/l). There is poor correlation between turbidity and triglyceride concentration.

Rheumatoid factors < 1200 IU/ml do not interfere.

A high-dose hook effect may occur at CRP concentrations > 50 mg/dl.

Measuring/reportable range

Roche/Hitachi 704/902

Measuring range: 0.3–24 mg/dl (0.003–0.24 g/l)*

When the CRP concentration of the sample is above the measuring range, manually dilute with 0.9% NaCl or CRP-free serum (e.g. 1 + 2). Multiply the result by the appropriate dilution factor (e.g. 3).

Roche/Hitachi 717/904/911/912/914/917/MODULAR

Measuring range: 0.3–24 mg/dl (0.003–0.24 g/l)*

Extended measuring range with rerun: 0.3–48 mg/dl (0.003–0.48 g/l)**

Roche/Hitachi 747

Measuring range: 0.3–20 mg/dl (0.003–0.20 g/l)*

Extended measuring range with rerun: 0.3–40 mg/dl (0.003–0.40 g/l)**

* The maximum reportable range is dependent on the highest standard concentration.

** The maximum extended reportable range is an approximate value dependent on the calculated value of the highest standard.

Expected values

IFCC/CRM 470***		Roche	
mg/dl	g/l	mg/dl	mg/l
< 0.5	< 0.005	< 0.5	< 5.0

***Reference values according to CRM 470 protein standardization¹¹

Elevated CRP range¹²

Concentrations of < 1 mg/dl (< 0.01 g/l) exclude many acute inflammatory diseases, but do not specifically exclude inflammatory processes.

Elevated concentrations of < 5 mg/dl (< 0.05 g/l) in acute disease occur in the presence of slight to moderate inflammatory processes. Values > 5 mg/dl (> 0.05 g/l) indicate high and extensive inflammatory activity.

It is important to monitor the CRP concentration during the acute phase of the illness.

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference range. For diagnostic purposes, CRP results should always be assessed in conjunction with the patient's medical history, clinical examinations and other findings.

Specific performance data

The data determined using a Roche/Hitachi system are given below. Results obtained in individual laboratories may differ.

Imprecision¹⁰

Reproducibility was determined using human samples and controls in an internal protocol (within run n = 21, between day n = 8). The following results were obtained.

Sample	Within run			Between day		
	Mean mg/dl	SD mg/dl	%CV	Mean mg/dl	SD mg/dl	%CV
Human serum 1	2.3	0.03	1.3	2.2	0.13	6.0
Human serum 2	5.4	0.05	1.0	4.2	0.12	2.9
Human serum 3	9.4	0.06	0.6	11.5	0.15	1.3

Analytical sensitivity (lower detection limit)¹⁰

0.3 mg/dl (0.003 g/l)

The lower detection limit represents the lowest measurable CRP concentration that can be distinguished from zero. It is calculated as three standard deviations of 21 replicates of the lowest standard.

Method comparison¹⁰

A comparison of the CRP determination using the Roche Tina-quant \square CRP assay on a Roche/Hitachi 911 analyzer (y) with a nephelometric method (x) gave the following correlation in human sera (mg/dl):

Passing/Bablok^{13,14} Linear regression
y = 0.179 + 1.051 x y = 0.292 + 1.016 x
r = 0.991 r = 0.991

SD (md 95) = 0.758 Sy.x = 0.517

Number of samples measured: 85

The sample concentrations were between 0.10 and 23.6 mg/dl.

References

- 1 Henry JB, ed. Clinical Diagnosis and Management by Laboratory Methods, Vol II. Philadelphia, Pa: WB Saunders Co, 1979.
- 2 Greiling H, Gressner AM, eds. Lehrbuch der Klinischen Chemie und Pathobiochemie, 3rd ed. Stuttgart/New York: Schattauer, 1995: 234–236.
- 3 Thomas L, Messenger M. Pathobiochemie und Labordiagnostik der Entzündung. Lab med 1993;17:179–194.
- 4 Young B, Gleeson M, Cripps AW. C-reactive protein: A critical review. Pathology 1991;23:118–124.
- 5 Tietz NW. Clinical Guide to Laboratory Tests, 3rd ed Philadelphia, Pa: WB Saunders, 1995.
- 6 Lizana J, Hellsing K. Clin Chem 1974;20:1181.
- 7 Tietz NW. Fundamentals of Clinical Chemistry, 2nd ed Philadelphia, Pa: WB Saunders, 1976:278–280.
- 8 Heidelberger M, Kendall FE. J Exp Med 1935;62:697.
- 9 Gunder WG, Narayanan S, Wisser H, Zawta B. List of Analytes; Pre-analytical Variables. Brochure in: Samples: From the Patient to the Laboratory. Darmstadt: GIT Verlag, 1996.
- 10 Data on file at Roche.
- 11 Konsensuswerte der Deutschen Gesellschaft für Laboratoriumsmedizin, der Deutschen Gesellschaft für Klinische Chemie und des Verbandes der Diagnostica-Industrie e.V. (VDGH). Clin Lab 1995;41:743–748.
- 12 Thomas L, ed. Labor und Diagnose, 4th ed Marburg: Die Medizinische Verlagsgesellschaft, 1992:781.
- 13 Passing H, Bablok W. A New Biometrical Procedure for Testing the Equality of Measurements from Two Different Analytical Methods. J Clin Chem Clin Biochem 1983;21:709–720.
- 14 Bablok W et al. A General Regression Procedure for Method Transformation. J Clin Chem Clin Biochem 1988;26:783–790.

Instrument settings

US users

Refer to application sheet for additional operating instructions.

Roche/Hitachi 914 customers

Refer to application sheet for parameters.

Roche/Hitachi 904, 911, 912, 917, MODULAR users

Read in the application parameters from the application diskette, barcode sheet or settings sheet as appropriate.

Roche/Hitachi 704

Temperature: 37°C

PROGRAM 2 CHEMISTRY PARAMETERS	
TEST	[CRP]
ASSAY CODE	[2(2 POINT)]-[15]-[32]
SAMPLE VOLUME	[20]
R1 VOLUME	[350]-[20]-[NO]
R2 VOLUME	[50]-[20]-[NO]
WAVELENGTH	[700]-[340]
CALIB. METHOD	[NONLINEAR]-[1]-[6]
STD. (1) CONC.-POS.	[]-[]
STD. (2) CONC.-POS.	[]-[]
STD. (3) CONC.-POS.	[]-[]
STD. (4) CONC.-POS.	[]-[]
STD. (5) CONC.-POS.	[]-[]
STD. (6) CONC.-POS.	[]-[]
UNIT	[]
SD LIMIT	[80]
DUPLICATE LIMIT	[200]
SENSITIVITY LIMIT	[2200]
ABS. LIMIT (INC/DEC)	[0]-[INCREASE]
PROZONE LIMIT	[32000]-[UPPER]
EXPECTED VALUES	[]-[]
INSTRUMENT FACTOR	[1.00]

— Data entered by the operator

Roche/Hitachi 717

Temperature: 37°C

PROGRAM 2 CHEMISTRY PARAMETERS	
TEST	[CRP]
ASSAY CODE	[2(2 POINT)]-[24]-[50]
SAMPLE VOLUME	[20]-[10]
R1 VOLUME	[350]-[100]-[NO]
R2 VOLUME	[50]-[50]-[NO]
WAVELENGTH	[700]-[340]
CALIB. METHOD	[NONLINEAR]-[1]-[6]
STD. (1) CONC.-POS.	[]-[]
STD. (2) CONC.-POS.	[]-[]
STD. (3) CONC.-POS.	[]-[]
STD. (4) CONC.-POS.	[]-[]
STD. (5) CONC.-POS.	[]-[]
STD. (6) CONC.-POS.	[]-[]
SD LIMIT	[80]
DUPLICATE LIMIT	[200]
SENSITIVITY LIMIT	[2200]
ABS. LIMIT (INC/DEC)	[0]-[INCREASE]
PROZONE LIMIT	[32000]-[UPPER]
EXPECTED VALUES	[]-[]
PANIC VALUES	[]-[]
INSTRUMENT FACTOR	[1.00]

— Data entered by the operator

Roche/Hitachi 747

Temperature: 37°C

PROGRAM 4.2 CHEMISTRY PARAMETERS		
TEST	[CRP]	
ASSAY CODE	[2(2POINT)]-[22]-[50]	
WAVELENGTH	[700 (SUB)]-[340 (MAIN)]	
	SERUM	URINE
SAMPLE VOLUME	[20]-[10]	[]-[]
EXPECTED VALUES	[]-[]	[]-[]
PANIC VALUES	[]-[]	[]-[]
ABS. LIMIT (INC/DEC)	[0]-[INCREASE]	[]-[]
PROZONE LIMIT	[32000]-[UPPER]	[]-[]
	R1	R2
R1/R2 VOLUME	[350]	[50]
R1/R2 DUMMY INTERVAL	[0]	[0]
DILUTION VOLUME	[0]	
CALIBRATION TYPE	[NONLINEAR]	
POINTS	[6]	
STD 1 CONC RACK POS	[]-[]-[]	
STD 2 CONC RACK POS	[]-[]-[]	
STD 3 CONC RACK POS	[]-[]-[]	
STD 4 CONC RACK POS	[]-[]-[]	
STD 5 CONC RACK POS	[]-[]-[]	
STD 6 CONC RACK POS	[]-[]-[]	
SD LIMIT	[120]	
DUPLICATE LIMIT	[200]	
SENSITIVITY LIMIT	[2500]	
STD 1 ABS. LEVEL	[]-[]	
INSTRUMENT FACTOR	[1.00]	

— Data entered by the operator

Roche/Hitachi 902


No.	<Chemistry>	
1	Test Name	CRP
2	Assay Code (Mthd)	2 Point End
3	Assay Code (2.Test)	0
4	Reaction Time	10
5	Assay Point 1	17
6	Assay Point 2	35
7	Assay Point 3	0
8	Assay Point 4	0
9	Wavelength. (SUB)	700
10	Wavelength (MAIN)	340
11	Sample Volume	20.0
12	R1 Volume	350
13	R1 Pos.
14	R1 Bottle Size	Small
15	R2 Volume	0
16	R2 Pos.	0
17	R2 Bottle Size	Small
18	R3 Volume	50
19	R3 Pos.
20	R3 Bottle Size	Small
21	Calib. Type (Type)	Logit-Log(4P)
22	Calib. Type (Wght)	0
23	Calib. Conc. 1	0.0
24	Calib. Pos. 1
25	Calib. Conc. 2
26	Calib. Pos. 2
27	Calib. Conc. 3
28	Calib. Pos. 3
29	Calib. Conc. 4
30	Calib. Pos. 4
31	Calib. Conc. 5
32	Calib. Pos. 5
33	Calib. Conc. 6
34	Calib. Pos. 6
35	S1 ABS.	0
36	K Factor	10000
37	K2 Factor	10000
38	K3 Factor	10000
39	K4 Factor	10000
40	K5 Factor	10000
41	A Factor	0
42	B Factor	0
43	C Factor	0
44	SD Limit	80
45	Duplicate Limit	200
46	Sens. Limit	2500
47	S1ABS. Limit (L)	-32000
48	S1ABS. Limit (H)	32000
49	ABS. Limit	0
50	ABS. Limit (D/I)	Increase
51	Prozone Limit	32000
52	Proz Limit (Upp/Lower)	Upper
53	Prozone (End Point)	35
54	Expect. Value (L)
55	Expect. Value (H)
56	Instr. Fact. (a)	1.0
57	Instr. Fact. (b)	0.0
58	Key Setting

.... Data entered by the operator


For detailed information, consult the operator manuals for Roche/Hitachi systems, the respective application sheets and the package inserts for the calibrators and control sera. Tina-quant, Preciset, Precinorm and Precipath are trademarks of a member of the Roche Group.

Intralipid is a trademark of KabiPharmacia, Inc.

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Roche/Hitachi 747

Temperatura: 37°C

PROGRAM 4.2 CHEMISTRY PARAMETERS		
TEST	[CRP]	
ASSAY CODE	[2(2POINT)]-[22]-[50]	
WAVELENGTH	[700 (SUB)]-[340 (MAIN)]	
	SERUM	URINE
SAMPLE VOLUME	[20]-[10]	[]-[]
EXPECTED VALUES	[]-[]	[]-[]
PANIC VALUES	[]-[]	[]-[]
ABS. LIMIT (INC/DEC)	[0]-[INCREASE]	[]-[]
PROZONE LIMIT	[32000]-[UPPER]	[]-[]
	R 1	R 2
R1/R2 VOLUME	[350]	[50]
R1/R2 DUMMY INTERVAL	[0]	[0]
DILUTION VOLUME	[0]	
CALIBRATION TYPE	[NONLINEAR]	
POINTS	[6]	
STD 1 CONC RACK POS	[]-[]-[]	
STD 2 CONC RACK POS	[]-[]-[]	
STD 3 CONC RACK POS	[]-[]-[]	
STD 4 CONC RACK POS	[]-[]-[]	
STD 5 CONC RACK POS	[]-[]-[]	
STD 6 CONC RACK POS	[]-[]-[]	
SD LIMIT	[120]	
DUPLICATE LIMIT	[200]	
SENSITIVITY LIMIT	[2500]	
STD 1 ABS. LEVEL	[]-[]	
INSTRUMENT FACTOR	[1.00]	

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©2001 Roche Diagnostics
Roche Diagnostics GmbH, D-68298 Mannheim, Alemanha