## Lactate

• Indicates analyzer(s) on which kit(s) can be used

Cat. no.	Bottle	Contents	704	717	736 737	747	902	904	911 912	914	917	MODULAR	
												P	D
1822837	1	Buffer, Enzymes 2 x 20 ml	•	•		•	•	•	•	•	•	•	
	2	Buffer, Enzymes 2 x 6 ml											

Some analyzers and kits shown may not be available in all countries. For additional system applications, contact your local Roche Diagnostics representative.

### Intended use

For the quantitative determination of L-lactate in plasma, cerebrospinal fluid or whole blood manually and on automated clinical chemistry analyzers.

Anaerobic glycolysis markedly increases blood lactate and causes some increase in pyruvate levels, especially with prolonged exercise. The common cause for increased blood lactate and pyruvate is anoxia resulting from such conditions as shock, pneumonia and congestive heart failure. Lactic acidosis may also occur in renal failure and leukemia. Thiamine deficiency and diabetic ketoacidosis are associated with increased levels of lactate and pyruvate.

Lactate levels in cerebrospinal fluid (CSF) are increased in bacterial meningitis. Increased CSF levels also occur in hypocapnia, hydrocephalus, brain abscesses, cerebral ischemia and any clinical condition associated with reduced oxygenation of the brain and/or increased intracranial pressure.

Lactate measurements that evaluate the acid-base status are used in the diagnosis and treatment of lactic acidosis (abnormally high acidity in the blood).

### Summary

In recent years, enzymatic methods for the determination of lactate have gained favor over colorimetric and titrimetric methods. Enzymatic methods are generally simple and provide greater specificity, accuracy, and reproducibility.


The first enzymatic method described for the determination of lactate was based on the transfer of hydrogen from lactate to potassium ferricyanide by lactate dehydrogenase. However, the procedure was cumbersome and did not receive wide acceptance.

Subsequent methods involved the UV measurement of the formation of NADH. In 1974, Gutmann and Wahlefeld<sup>1</sup> described a lactate procedure that measures the NADH formed by the oxidation of lactate catalyzed by LD, using hydrazine as a trapping agent for pyruvate. A method described by Noll<sup>2</sup> is also based on the catalytic action of LD but includes ALT in the reaction mixture to more rapidly remove the pyruvate formed from the conversion of lactate.


The method presented here uses an enzymatic reaction to convert lactate to pyruvate. The hydrogen peroxide produced by this reaction is then used in an enzymatic reaction to generate a colored dye.<sup>3,4</sup> This method offers longer reagent stability than the previous UV enzymatic methods.

### Test principle

L-lactate is oxidized to pyruvate by the specific enzyme lactate oxidase (LOD).

Peroxidase (POD) is used to generate a colored dye using the hydrogen peroxide generated in the first reaction.<sup>3,4</sup>

The intensity of the color formed is proportional to the L-lactate concentration.

### Working solution concentration

**R1** Buffer, Enzymes  
Hydrogen donor; Ascorbate oxidase (cucumber): 30 U/ml; Buffers; Preservatives

**R2** Buffer, Enzymes  
4-aminoantipyrene: 1 mg/ml; LOD (microorganism): 15 U/ml; Peroxidase (horseradish): 24 U/ml; Buffers; Preservatives

### Precautions and warnings

For in vitro diagnostic use.  
Exercise the normal precautions required for handling all laboratory reagents.

### Reagent handling

R1: Ready for use.  
R2: Ready for use.

### Storage and stability

Unopened kit components: up to the expiration date at 2–8°C  
R1: 90 days opened and refrigerated on the analyzer  
R2: 90 days opened and refrigerated on the analyzer

### Specimen collection and preparation

Do not use serum. Use plasma from blood collected by standard venipuncture technique in fluoride-oxalate tubes (2.5 mg sodium fluoride and 2.0 mg potassium oxalate/ml blood). Centrifuge within 15 minutes of drawing specimen. Cerebrospinal fluid (CSF) may be used as obtained. Use whole blood collected with sodium fluoride/potassium oxalate anticoagulant.

#### NOTES:

- The lactate level increases rapidly with physical exercise. The time required for return to normal lactate values depends on the physical fitness of the subject. Thirty minutes at rest is usually sufficient for this purpose.
- Blood samples should be drawn from a stasis-free vein. However, minimal hemostasis (less than 30 seconds) will not affect lactate levels. Avoid the use of a tourniquet, if possible.<sup>5</sup>
- Glycolysis in blood samples can rapidly increase lactate levels. Cells contribute to the glycolysis and their quick removal is essential for accurate lactate analysis.<sup>6</sup> Heparinized plasma is acceptable, but precautions must be taken to retard glycolysis by keeping the whole blood on ice and then separating the plasma from the cells within 15 minutes of collection.<sup>8</sup>

### Whole blood must be deproteinized by the following procedure prior to assay:

Pipette into polypropylene microfuge tube:  
Trichloroacetic acid (TCA) 200 µl  
(10% w/v)  
Whole blood sample 200 µl

Stopper tightly, vortex, and centrifuge at 1500 RCF for 10 minutes. Remove stopper.

Pipette into centrifuged microfuge tube:  
0.9% saline 100 µl

Stopper tightly and invert tube gently several times. Transfer supernatant into a clean sample cup for analysis. Supernatant should be clear to slightly cloudy and colorless. Some particulate matter may be observed; this should settle to the bottom of the sample cup within several minutes. Multiply the whole blood supernatant result by 2.5 in order to correct for dilution.

**Stability:** Plasma(separated): 2 hours at 20–25 °C or 2 days at 2–8 °C.<sup>7</sup>  
CSF: 3 hours at 20–25 °C,  
24 hours at 2–8 °C or  
1 month at –20 °C.

### Testing procedure

#### Materials provided

- Working solutions as described above

#### Additional materials required

- Calibrators and controls as indicated below
- 0.9% NaCl

#### Assay

Refer to the appropriate operator's manual and/or the Instrument Settings section of this insert for analyzer specific assay instructions. The performance of applications not validated by Roche Diagnostics is not warranted and must be defined by the user.

#### Calibration

Standardization: the method was standardized against a gravimetrically prepared calibrator.

- S1: 0.9% NaCl or  
Precical Ethyl Alcohol/Ammonia/Lactate Calibrator Set, Bottle 1  
S2: Precical Calibrator Serum or  
Precical Ethyl Alcohol/Ammonia/Lactate Calibrator Set, Bottle 2 or  
C.f.a.s. (Calibrator for automated systems)

# Lactate

## Calibration frequency

2 point calibration is recommended

- after bottle change
  - after lot change
  - as required following quality control procedures
- Calibration verification: Not necessary.

## Quality control

For quality control use Precitrol-N, Precitrol-A, Precinorm U, Precipath U, Precitrol Ethyl Alcohol/Ammonia/Lactate Controls or other suitable control material. The control intervals and limits must be adapted to the individual laboratory and country-specific requirements. Values obtained should fall within established limits. Each laboratory should establish corrective measures to be taken if values fall outside the limits.

## Calculation

Roche/Hitachi systems automatically calculate the lactate concentration of each sample.

Conversion factor<sup>5</sup>: mg/dl x 0.111 = mmol/l

## Limitations - interference

Criterion: recovery ± 10% of initial values

Plasma

Icterus: No significant interference from unconjugated bilirubin up to an I index of 60 and conjugated bilirubin up to an I index of 28 (approximate unconjugated and conjugated bilirubin concentrations: 60 mg/dl and 28 mg/dl).

Hemolysis: No significant interference from hemoglobin up to an H index of 1000 (approximate hemoglobin concentration: 1000 mg/dl).

Lipemia (Intralipid): No significant interference from lipemia up to an L index of 1000 (approximate triglycerides concentration: 2000 mg/dl).

There is poor correlation between turbidity and triglycerides concentration.

Thirty-five commonly used pharmaceuticals were tested in vitro. Dopamine (10 mg/l), Levodopa (20 mg/l) and Methyldopa (20 mg/l) significantly reduced the lactate results. However, Dopamine at 1 mg/l, Levodopa at 4 mg/l and Methyldopa at 2 mg/l do not significantly affect lactate results. No interference with the assay was found with any of the other drugs tested.

## Measuring/reportable range

2.0 - 140 mg/dl (0.2 - 15.5 mmol/l)

Specimen dilution

Determine samples with lactate concentrations > 140 mg/dl via the rerun function. On instruments without rerun function, manually dilute samples with 0.9% NaCl (e.g. 1+1). Multiply the result by the appropriate dilution factor (e.g. 2).

## Expected values

Plasma:	4.5 - 19.8 mg/dl	(0.5 - 2.2 mmol/l)	venous <sup>5</sup>
CSF:	10 - 60 mg/dl	(1.1 - 6.7 mmol/l)	neonate <sup>5</sup>
	10 - 40 mg/dl	(1.1 - 4.4 mmol/l)	3 -10 days old
	10 - 25 mg/dl	(1.1 - 2.8 mmol/l)	> 10 days old
	10 - 22 mg/dl	(1.1 - 2.4 mmol/l)	adult
Whole Blood:	8.1 - 15.3 mg/dl	(0.9 - 1.7 mmol/l)	venous <sup>5</sup>
	< 11.3 mg/dl	(<1.3 mmol/l)	arterial

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference range.

For diagnostic purposes, lactate results should always be assessed in conjunction with the patient's medical history, clinical examinations and other findings.

## Specific performance data<sup>8</sup>

The data determined using a Roche/Hitachi system are given below. Results obtained in individual laboratories may differ.

## Imprecision

Reproducibility was determined for plasma using human samples and controls in an internal protocol (within run n=21; between day n=63). The following results were obtained.

Sample	Within run			Between day		
	Mean	SD	%CV	Mean	SD	%CV
	mg/dl	mg/dl		mg/dl	mg/dl	
Human plasma	17.6	0.07	0.4	17.8	0.18	1.0
Control 1	14.2	0.09	0.6	14.3	0.19	1.3
Control 2	38.7	0.10	0.3	39.0	0.44	1.1

Reproducibility was determined for CSF using CSF controls in an internal protocol (within run n=21; between day n=20). The following results were obtained.

Sample	Within run			Between day		
	Mean	SD	%CV	Mean	SD	%CV
	mg/dl	mg/dl		mg/dl	mg/dl	
Control 1	17.2	0.07	0.4	17.3	0.12	0.7
Control 2	62.4	0.36	0.6	62.8	0.46	0.7

## Analytical sensitivity (lower detection limit)

Detection limit: 2 mg/dl (0.2 mmol/l)

The lower detection limit represents the lowest measurable lactate concentration that can be distinguished from zero. It is calculated as three standard deviations of 21 replicates of the lowest standard.

## Method comparison

A comparison of the lactate determination using the Roche Diagnostics lactate assay (y) with another lactate determination (x) gave the following correlation (mg/dl):

### Plasma:

Passing/Bablok<sup>9,10</sup>

$y = -0.11 + 1.015x$

$r = 1.00$

SD (md 95) = 1.86

Number of samples measured: 37

The sample concentrations were between 7.0 and 133.9 mg/dl.

Linear regression

$y = 0.01 + 1.015x$

$r = 1.00$

$Sy.x = 1.20$

### CSF:

Passing/Bablok<sup>9,10</sup>

$y = 0.30 + 1.012x$

$r = 1.00$

SD (md 95) = 1.29

Number of samples measured: 42

The sample concentrations were between 11.8 and 129.1 mg/dl.

Linear regression

$y = 0.32 + 1.011x$

$r = 1.00$

$Sy.x = 0.92$

### Whole Blood:

Passing/Bablok<sup>9,10</sup>

$y = 0.04 + 1.036x$

$r = 0.997$

SD (md 95) = 2.9

Number of samples measured: 29

The sample concentrations were between 5.5 and 90.5 mg/dl.

Linear regression

$y = 0.13 + 1.040x$

$r = 0.997$

$Sy.x = 2.26$

## References

- 1 Gutmann I, Wahlefeld A. *Methods of Enzymatic Analysis*. 2nd ed. Bergmeyer HU (ed), New York, NY: Academic Press Inc: 1974;1464.
- 2 Noll F. *Methods of Enzymatic Analysis*. 2nd ed. Bergmeyer HU ed. New York, NY: Academic Press Inc; 1974;1475.
- 3 Trinder P. Determination of glucose in blood using glucose oxidase with an alternative oxygen acceptor. *Ann. Clin. Biochem.* 1969;6:24.
- 4 Barhan D, Trinder P. An improved colour reagent for the determination of blood glucose by the oxidase system. *Analyst* 97.1972:142.
- 5 Tietz NW. *Clinical Guide to Laboratory Tests*. 3rd ed. Philadelphia, PA: WB Saunders Co; 1995: 382-383.
- 6 Burtis CA, Ashwood ER, eds. *Tietz Textbook of Clinical Chemistry*. 2nd ed. Philadelphia, PA: W.B. Saunders Co; 1994: 976.
- 7 Westgard JO, Lahmeyer BL, Birnbaum ML. *Clin Chem.* 1972;18:1334-1338.
- 8 Data on file at Roche Diagnostics.
- 9 Passing H, Bablok W. A new biometrical procedure for testing the equality of measurements from two different analytical methods. *J Clin Chem Clin Biochem.* 1983;21:709-720.
- 10 Bablok W et al. A general regression procedure for method transformation. *J Clin Chem Biochem.* 1988;26:783-790.

# Lactate

## Instrument settings

US users: Refer to application sheet for additional operating information.  
 Roche/Hitachi 914 customers:  
 Refer to application sheet for parameters.  
 Roche/Hitachi 904, 911, 912, 917, MODULAR users:  
 Read in the application parameters from the application diskette or barcode sheet, as appropriate.

### Roche/Hitachi 704

Temperature: 37°C

PROGRAM 2 CHEMISTRY PARAMETERS	
TEST	[LAC]
ASSAY CODE	[2 POINT] - [15] - [23]
SAMPLE VOLUME	[4]
R1 VOLUME	[350] - [20] - [NO]
R2 VOLUME	[ 70] - [20] - [NO]
WAVELENGTH	[700] - [660]
CALIB. METHOD	[LINEAR] - [0] - [0]
STD. (1) CONC.-POS.	[ .... ] - [ .... ]
STD. (2) CONC.-POS.	[ .... ] - [ .... ]
STD. (3) CONC.-POS.	[0] - [0]
STD. (4) CONC.-POS.	[0] - [0]
STD. (5) CONC.-POS.	[0] - [0]
STD. (6) CONC.-POS.	[0] - [0]
UNIT	[ .... ]
SD LIMIT	[0.1]
DUPLICATE LIMIT	[ 25]
SENSITIVITY LIMIT	[100]
ABS. LIMIT (INC/DEC)	[32000] - [INCREASE]
PROZONE LIMIT	[0] - [LOWER]
EXPECTED VALUE	[ .... ] - [ .... ]
INSTRUMENT FACTOR	[1.00]

.... Data entered by the operator

### Roche/Hitachi 747

Temperature: 37°C

PROGRAM 4.2 CHEMISTRY PARAMETERS		
TEST	..... (LAC)	
ASSAY CODE	2 POINT - 22 - 34	
WAVELENGTH	700 (SUB) - 660 (MAIN)	
	<b>SERUM</b>	<b>URINE</b>
SAMPLE VOLUME (µl)	3 - 2	[ 3 ] - [ 2 ]
EXPECTED VALUE	..... - .....	[.....] - [.....]
PANIC VALUE	..... - .....	[.....] - [.....]
ABS. LIMIT (INC/DEC)	32000 - INC	[32000] - [INC]
PROZONE LIMIT	0 - LOWER	[0]-[LOWER]
	<b>R1</b>	<b>R2</b>
R1/R2 VOLUME (µl)	250	50
R1/R2 DUMMY INTERVAL	0	0
DILUTION VOLUME (µl)	0	
CALIB. METHOD	LINEAR	
POINTS	0	
STD 1 CONC RACK POS (NaCl)	..... - ..... - .....	
STD 2 CONC RACK POS	..... - ..... - .....	
STD 3 CONC RACK POS	0 - 0 - 0	
STD 4 CONC RACK POS	0 - 0 - 0	
STD 5 CONC RACK POS	0 - 0 - 0	
STD 6 CONC RACK POS	0 - 0 - 0	
SD LIMIT	0.1	
DUPLICATE LIMIT	25	
SENSITIVITY LIMIT	100	
STD 1 ABS. LEVEL	-100 - 100	
INSTRUMENT FACTOR	1.0	

.... Data entered by the operator

### Roche/Hitachi 717

Temperature: 37°C

PROGRAM 2 CHEMISTRY PARAMETERS	
TEST	[LAC]
ASSAY CODE	[2 POINT] - [24] - [36]
SAMPLE VOLUME	[3] - [2]
R1 VOLUME	[250] - [20] - [NO]
R2 VOLUME	[50] - [20] - [NO]
WAVELENGTH	[700] - [660]
CALIB. METHOD	[LINEAR] - [0] - [0]
STD. (1) CONC.-POS.	[... ] - [... ]
STD. (2) CONC.-POS.	[... ] - [... ]
STD. (3) CONC.-POS.	[0] - [0]
STD. (4) CONC.-POS.	[0] - [0]
STD. (5) CONC.-POS.	[0] - [0]
STD. (6) CONC.-POS.	[0] - [0]
SD LIMIT	[0.1]
DUPLICATE LIMIT	[25]
SENSITIVITY LIMIT	[100]
ABS. LIMIT (INC/DEC)	[32000] - [INCREASE]
PROZONE LIMIT	[0] - [LOWER]
EXPECTED VALUE	[ .... ] - [ .... ]
PANIC VALUE	[ .... ] - [ .... ]
INSTRUMENT FACTOR	[1.00]

.... Data entered by the operator

# Lactate

## Roche/Hitachi 902

No.	<Chemistry>	
1	Test Name	LAC
2	Assay Code (Mthd)	2 Point End
3	Assay Code (2. Test)	0
4	Reaction Time	5
5	Assay Point 1	4
6	Assay Point 2	15
7	Assay Point 3	0
8	Assay Point 4	0
9	Wave Leng. (SUB)	700
10	Wave Leng. (MAIN)	660
11	Sample Volume	3.0
12	R1 Volume	250
13	R1 Pos.	.....
14	R1 Bottle Size	Small
15	R2 Volume	50
16	R2 Pos.	.....
17	R2 Bottle Size	Small
18	R3 Volume	0
19	R3 Pos.	0
20	R3 Bottle Size	Small
21	Calib. Type (Type)	Linear
22	Calib. Type (Wght)	0
23	Calib. Conc. 1	.....
24	Calib. Pos. 1	.....
25	Calib. Conc. 2	.....
26	Calib. Pos. 2	.....
27	Calib. Conc. 3	0
28	Calib. Pos. 3	0
29	Calib. Conc. 4	0
30	Calib. Pos. 4	0
31	Calib. Conc. 5	0
32	Calib. Pos. 5	0
33	Calib. Conc. 6	0
34	Calib. Pos. 6	0
35	S1 ABS.	0
36	K Factor	10000
37	K2 Factor	10000
38	K3 Factor	10000
39	K4 Factor	10000
40	K5 Factor	10000
41	A Factor	0
42	B Factor	0
43	C Factor	0
44	SD Limit	0.1
45	Duplicate Limit	25
46	Sens. Limit	100
47	S1ABS. Limit (L)	-100
48	S1ABS. Limit (H)	100
49	ABS. Limit	32000
50	ABS. Limit (D/I)	Increase
51	Prozone Limit	0
52	Proz Limit (Upp/Low)	Lower
53	Proz Limit (End Point)	35
54	Expect. Value (L)	.....
55	Expect. Value (H)	.....
56	Instr. Fact. (a)	1.0
57	Instr. Fact. (b)	0.0
58	Key Setting	.....

.... Data entered by the operator

For detailed information, consult the operator manuals for Roche/Hitachi systems, the respective application sheets and the package inserts for the calibrators and control sera.

©1998 Roche Diagnostics.  
050214304 04  
November 1998

Roche Diagnostics GmbH, Mannheim, Germany  
Roche Diagnostics Corporation, Indianapolis, IN, USA  
US Customer Technical Support 1-800-428-2336

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## Bibliographie

- 1 Gutmann I, Wahlefeld A. *Methods of Enzymatic Analysis*. 2nd ed. Bergmeyer HU (ed), New York, NY: Academic Press Inc; 1974;1464.
- 2 Noll F. *Methods of Enzymatic Analysis*. 2nd ed. Bergmeyer HU ed. New York, NY: Academic Press Inc; 1974;1475.
- 3 Trinder P, Determination of glucose in blood using glucose oxidase with an alternative oxygen acceptor. *Ann. Clin. Biochem.* 1969;6:24.
- 4 Barhan, D, Trinder, P. An improved colour reagent for the determination of blood glucose by the oxidase system. *Analyst* 97.1972:142.
- 5 Tietz NW. *Clinical Guide to Laboratory Tests*. 3rd ed. Philadelphia, PA: WB Saunders Co; 1995: 382-383.
- 6 Burtis CA, Ashwood ER, eds. *Tietz Textbook of Clinical Chemistry*. 2nd ed. Philadelphia, PA: W.B. Saunders Co; 1994: 976.
- 7 Westgard JO, Lahmeyer BL, Birnbaum ML. *Clin Chem.* 1972;18:1334-1338.
- 8 Dokumentation Roche Diagnostics.
- 9 Passing H, Bablok W. A new biometrical procedure for testing the equality of measurements from two different analytical methods. *J Clin Chem Clin Biochem.* 1983;21:709-720.
- 10 Bablok W et al. A general regression procedure for method transformation. *J Clin Chem Biochem.* 1988;26:783-790.

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- 1 Gutmann I, Wahlefeld A. *Methods of Enzymatic Analysis*. 2nd ed. Bergmeyer HU (ed), New York, NY: Academic Press Inc; 1974:1464.
- 2 Noll F. *Methods of Enzymatic Analysis*. 2nd ed. Bergmeyer HU ed. New York, NY: Academic Press Inc; 1974:1475.
- 3 Trinder P. Determination of glucose in blood using glucose oxidase with an alternative oxygen acceptor. *Ann. Clin. Biochem.* 1969;6:24.
- 4 Barhan D, Trinder P. An improved colour reagent for the determination of blood glucose by the oxidase system. *Analyst* 97.1972:142.
- 5 Tietz NW. *Clinical Guide to Laboratory Tests*. 3rd ed. Philadelphia, PA: WB Saunders Co; 1995: 382-383.
- 6 Burtis CA, Ashwood ER, eds. *Tietz Textbook of Clinical Chemistry*. 2nd ed. Philadelphia, PA: W.B. Saunders Co; 1994: 976.
- 7 Westgard JO, Lahmeyer BL, Birnbaum ML. *Clin Chem.* 1972;18:1334-1338.
- 8 Data on file at Roche Diagnostics.
- 9 Passing H, Bablok W. A new biometrical procedure for testing the equality of measurements from two different analytical methods. *J Clin Chem Clin Biochem.* 1983;21:709-720.
- 10 Bablok W et al. A general regression procedure for method transformation. *J Clin Chem Biochem.* 1988;26:783-790.

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## Riferimenti

- 1 Gutmann I, Wahlefeld A. *Methods of Enzymatic Analysis*. 2nd ed. Bergmeyer HU (ed), New York, NY: Academic Press Inc; 1974:1464.
- 2 Noll F. *Methods of Enzymatic Analysis*. 2nd ed. Bergmeyer HU ed. New York, NY: Academic Press Inc; 1974:1475.
- 3 Trinder P. Determination of glucose in blood using glucose oxidase with an alternative oxygen acceptor. *Ann. Clin. Biochem.* 1969;6:24.
- 4 Barhan D, Trinder P. An improved colour reagent for the determination of blood glucose by the oxidase system. *Analyst* 97.1972:142.
- 5 Tietz NW. *Clinical Guide to Laboratory Tests*. 3rd ed. Philadelphia, PA: WB Saunders Co; 1995: 382-383.
- 6 Burtis CA, Ashwood ER, eds. *Tietz Textbook of Clinical Chemistry*. 2nd ed. Philadelphia, PA: W.B. Saunders Co; 1994: 976.
- 7 Westgard JO, Lahmeyer BL, Birnbaum ML. *Clin Chem.* 1972;18:1334-1338.
- 8 Documentazione da Roche Diagnostics.
- 9 Passing H, Bablok W. A new biometrical procedure for testing the equality of measurements from two different analytical methods. *J Clin Chem Clin Biochem.* 1983;21:709-720.
- 10 Bablok W et al. A general regression procedure for method transformation. *J Clin Chem Biochem.* 1988;26:783-790.

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Roche Diagnostics GmbH, Germania