LDL-Cholesterol, no pretreatment

● Indicates Roche/Hitachi analyzer(s) on which kit(s) can be used

Cat. no.	Bottle	Contents	704	717	736 737	747	902	904	911 912	914	917	MODULAR	
												P	D
1985604	1	α -Cyclodextrin/buffer, 6 x 14 ml	●				●	●	●				
	2	Buffer/enzymes/4-aminoantipyrine, 3 x 10 ml											
1985612	1	α -Cyclodextrin/buffer, 6 x 36 ml		●				●	●	●			
	2	Buffer/enzymes/4-aminoantipyrine, 3 x 25 ml											
1985639	1	α -Cyclodextrin/buffer, 6 x 54 ml									●	●	
	2	Buffer/enzymes/4-aminoantipyrine, 6 x 20 ml											
1875892 1875914	1	α -Cyclodextrin/buffer, 6 x 267 ml										●	●
	2	Buffer/enzymes/4-aminoantipyrine, 6 x 103 ml											
1985647 1985655	1	α -Cyclodextrin/buffer, 5 x 252 ml			●	●							
	2	Buffer/enzymes/4-aminoantipyrine, 3 x 110 ml											

Some analyzers and kits shown may not be available in all countries. For additional system applications, contact your local Roche representative.

Intended use

Homogeneous enzymatic assay for the direct quantitative determination of LDL-cholesterol in human serum and plasma on automated clinical chemistry analyzers.

Summary

Low Density Lipoproteins (LDL) play a key role in causing and influencing the progression of atherosclerosis and coronary sclerosis in particular. The LDLs are derived from VLDLs (Very Low Density Lipoproteins) rich in triglycerides by the action of various lipolytic enzymes and are synthesized in the liver. The elimination of LDL from plasma takes place mainly by liver parenchymal cells via specific LDL receptors. Elevated LDL concentrations in blood and an increase in their residence time coupled with an increase in the biological modification rate results in the destruction of the endothelial function and a higher LDL-cholesterol uptake in the monocyte/macrophage system as well as by smooth muscle cells in vessel walls. The majority of cholesterol stored in atherosclerotic plaques originates from LDL.

The LDL-cholesterol value is the most powerful clinical predictor among all of the single parameters with respect to coronary atherosclerosis. Therefore, therapies focusing on lipid reduction primarily target the reduction of LDL-cholesterol which is then expressed in an improvement of the endothelial function, prevention of atherosclerosis and reducing its progression as well as preventing plaque rupture.

Various methods are available for the determination of LDL-cholesterol such as ultracentrifugation as the reference method, lipoprotein electrophoresis and precipitation methods. In the precipitation methods apolipoprotein-B-containing LDL-cholesterol is, for example, precipitated using either polyvinyl sulfate, dextran sulfate or polycyclic anions. The LDL-cholesterol content is usually calculated from the difference between total cholesterol and cholesterol in the remainder (VLDL- and HDL-cholesterol) in the supernate after precipitation with polyvinyl sulfate and dextran sulfate. Lipid Research Clinics recommend a combination of ultracentrifugation and precipitation methods using polyanions in the presence of divalent cations. The precipitation methods are however time-consuming, cannot be automated and are susceptible to interference by hyperlipidemic serum, particularly at high concentrations of free fatty acids. A more recent method is based on the determination of LDL-cholesterol after the sample is subjected to immunoabsorption and centrifugation.

The calculation of the LDL-cholesterol concentration according to Friedewald's formula is commonly practised. The formula is based on 2 cholesterol determinations, 1 triglyceride determination as well as precipitation of the HDL particles and presumes that a direct relationship exists between VLDL-cholesterol and triglycerides in fasting blood samples. Even in the presence of small amounts of chylomicrons or abnormal lipoproteins, the formula gives rise to falsely low LDL-cholesterol values. For this reason a great need exists for a simple and reliable method for the determination of LDL-cholesterol without any preparatory steps or calculation.


This automated method for the direct determination of LDL-cholesterol takes advantage of the selective micellar solubilization of LDL-cholesterol by a nonionic detergent and the interaction of a sugar compound and lipoproteins (VLDL and chylomicrons). When a detergent is included in the enzymatic method for cholesterol determination (cholesterol esterase cholesterol oxidase coupling reaction), the relative reactivities of cholesterol in the lipoprotein fractions increase in this order: HDL < chylomicrons < VLDL < LDL. In the presence of Mg⁺⁺, a sugar compound markedly reduces the enzymatic reaction of the cholesterol measurement in VLDL and chylomicrons. The combination of a sugar compound with detergent enables the selective determination of LDL-cholesterol in serum.¹⁻⁸

Nonfasting sample results are slightly lower than fasting results. Comparable non-fasting results were observed with the beta quantification method.⁹⁻¹¹ This direct assay meets the 1995 NECP goals of < 4% Total CV, bias ≤ 4% versus reference method, and ≤ 12% total analytical error.¹²





Test principle

Homogeneous enzymatic colorimetric assay

- Sample and addition of R1

- Addition of R2 and start of reaction:

Working solution concentration

R1 α -Cyclodextrin/buffer

MOPS (3-morpholinopropane sulfonic acid buffer): 50 mmol/l, pH 7.0; α -cyclodextrin sulfate: 0.1 g/l; dextran sulfate: 0.7 g/l; Mg₂SO₄·7H₂O: 7 g/l; HSDA^{**}: 0.3 g/l; AOD (Eupenicillium spec., recombinant): 3 U/l; peroxidase (horseradish): 10 U/l; preservative

R2 Buffer/enzymes/4-aminoantipyrine

MOPS: 50 mmol/l, pH 7.0; 4-aminoantipyrine: 0.51 g/l; cholesterol esterase (Pseudomonas spec.) > 0.7 U/l; cholesterol oxidase (Pseudomonas spec.) ≥ 0.3 U/l; peroxidase (horseradish) ≥ 20 U/l; detergent; preservative

* AA = 4-aminoantipyrine

** HSDA = Sodium N-(2-hydroxy-3-sulfopropyl)-3,5-dimethoxyaniline
AOD = ascorbate oxidase

Precautions and warnings

For in vitro diagnostic use.

Exercise the normal precautions required for handling all laboratory reagents.

Reagent handling

R1: Ready for use

R2: Ready for use

Storage and stability

Unopened kit components: Up to the expiration date at 2–8°C

R1: 28 days opened and refrigerated on the analyzer

R2: 28 days opened and refrigerated on the analyzer

LDL-C Plus

Specimen collection and preparation

Collect serum using standard sampling tubes.

Li- or Na-heparin plasma.

Stability¹¹: 7 days at 2–8°C

30 days at -70°C

Fasting and non-fasting samples can be used.¹⁰ EDTA plasma causes decreased values.

Centrifuge samples containing precipitate before performing the assay.

Testing procedure

Materials provided

- Working solutions as described above

Additional materials required

- Calibrators and controls as indicated below
- 0.9% NaCl

Assay

Refer to the appropriate operator's manual and/or the Instrument Settings section of this insert for analyzer-specific assay instructions. The performance of applications not validated by Roche is not warranted and must be defined by the user.

Calibration

Standardization: The method is traceable to the beta quantification method as defined in the recommendations in the LDL Cholesterol Method Certification Protocol for Manufacturers.¹²

S1: 0.9% NaCl

S2: C.f.a.s. (Calibrator for automated systems) HDL/LDL-C plus

Calibration frequency

Two point calibration is recommended

- after reagent lot change
- as required following quality control procedures

Calibration verification: Not necessary.

Quality control

For quality control, use Precinorm L, Precipath HDL/LDL-C, or other suitable control material.

The control intervals and limits must be adapted to the individual laboratory and country-specific requirements. Values obtained should fall within established limits. Each laboratory should establish corrective measures to be taken if values fall outside the limits.

Calculation

Roche/Hitachi systems automatically calculate the LDL-cholesterol concentration of each sample.

Conversion factor: $\text{mg/dl} \times 0.0259 = \text{mmol/l}$
 $\text{mmol/l} \times 38.66 = \text{mg/dl}$

Limitations – interference¹¹

(US users refer to application sheet for special wash instructions.)

To avoid probe and cuvette carry-over to Tina-quant \square Albumin, Lp(a), Lipase, and Magnesium, use the Evasion software with SMS/Acid Wash, Multiclean, or NaOH D in the LDL-cholesterol assay.

Criterion: Recovery within $\pm 10\%$ of initial values.

Icterus: No significant interference up to an I index of 60 (approximate conjugated and unconjugated bilirubin concentration: 60 mg/dl).

Hemolysis: No significant interference up to an H index of 1000 (approximate hemoglobin concentration: 1000 mg/dl).

Lipemia (Intralipid): No significant interference up to an L index of 200. Dilute samples with L index > 200 with 0.9% NaCl (e.g. 1 + 3).

No significant interference from native triglycerides up to 1200 mg/dl. Dilute samples with triglycerides > 1200 with 0.9% NaCl (e.g. 1 + 3). Multiply the result by the appropriate dilution factor (e.g. 4).

Rheumatoid factors < 200 IU/ml do not interfere. Ascorbic acid up to 50 mg/dl and oleic acid up to 4.0 mmol/l do not interfere.

No significant interference from HDL, VLDL, or chylomicrons.

In rare cases, elevated immunoglobulin concentrations can lead to falsely elevated LDL-cholesterol results.

Abnormal liver function does affect lipid metabolism; consequently HDL and LDL results are of limited diagnostic value. In some patients with abnormal liver function, the LDL-C plus result is significantly negatively biased versus beta quantification results.

Measuring/reportable range

3–550 mg/dl (0.03–5.5 g/l or 0.077–14.2 mmol/l)

Determine samples with LDL-cholesterol concentrations > 550 mg/dl via the rerun function. On instruments without rerun function, manually dilute samples with 0.9% NaCl (e.g. 1 + 1). Multiply the result by the appropriate dilution factor (e.g. 2).

Expected values^{8,13}

Levels in terms of risk for coronary heart disease:

Adult levels:

Recommended (desirable) < 130 mg/dl (<3.37 mmol/l)

Moderate risk: 130–159 mg/dl (3.37–4.12 mmol/l)

High risk: ≥ 160 mg/dl (≥ 4.14 mmol/l)

Recommended values according to the GRIPS study¹⁴

mg/dl	mmol/l	
145	3.8	For patients with manifested coronary heart disease
170	4.4	For patients having one or more risk factors
200	5.2	For persons exhibiting no risk factors and without manifested coronary heart disease

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference range. For diagnostic purposes, LDL-cholesterol results should always be assessed in conjunction with the patient's medical history, clinical examinations and other findings.

Specific performance data

The data determined using a Roche/Hitachi system are given below. Results obtained in individual laboratories may differ.

Imprecision¹¹

Reproducibility was determined using human samples and controls in an internal protocol: n = 21. The following results were obtained.

Sample	Within run			Between day		
	Mean mg/dl	SD mg/dl	%CV	Mean mg/dl	SD mg/dl	%CV
Human serum	108.2	0.83	0.8	108.4	2.93	2.7
Precinorm L	94.8	0.75	0.8	94.8	2.60	2.7
Precipath HDL/LDL-C	211.1	2.19	1.0	210.3	3.93	1.9

Analytical sensitivity (lower detection limit)¹¹

Detection limit: 3 mg/dl (0.03 g/l or 0.078 mmol/l)

The lower detection limit represents the lowest measurable LDL-cholesterol concentration that can be distinguished from zero. It is calculated as the concentration lying three standard deviations above that of the lowest standard (zero standard + 3SD, within run precision, n = 21).

Method comparison¹¹

- A. A comparison of the LDL-cholesterol determination using the Roche LDL-cholesterol plus assay (y) with another LDL-cholesterol determination (x) using a similar direct method on a Roche/Hitachi 911 analyzer gave the following correlation (mg/dl):
- | | |
|---------------------------------|------------------------|
| Passing Bablok ^{15,16} | Linear regression |
| $y = -12.46 + 1.041 x$ | $y = -13.50 + 1.051 x$ |
| $r = 0.978$ | $r = 0.980$ |
| $SD(\text{md } 95) = 12.76$ | $Sy.x = 6.62$ |
- Number of samples measured: 138
The sample concentrations were between 32.5 and 284.4 mg/dl.
- B. A comparison of the LDL-cholesterol determination using the Roche LDL-cholesterol plus assay on a Roche/Hitachi 911 analyzer (y) with the beta quantification method (ultra-centrifugation and precipitation with phosphotungstic acid/magnesium chloride) LDL-cholesterol determination (x) gave the following correlation (mg/dl):
- | | |
|---------------------------------|------------------------|
| Passing Bablok ^{15,16} | Linear regression |
| $y = -5.597 + 1.018 x$ | $y = -2.590 + 1.002 x$ |
| $r = 0.985$ | $r = 0.985$ |
| $SD(\text{md } 95) = 9.13$ | $Sy.x = 3.95$ |
- Number of samples measured: 40
The sample concentrations were between 69 and 208 mg/dl.
- C. A comparison on a Roche/Hitachi 911 analyzer of the LDL-cholesterol determination using the Roche LDL-cholesterol plus assay (y) with the PVS LDL-cholesterol determination (x) gave the following correlation (mg/dl):
- | | |
|---------------------------------|---------------------|
| Passing Bablok ^{15,16} | Linear regression |
| $y = 1.77 + 0.97 x$ | $y = 5.86 + 0.94 x$ |
| $r = 0.983$ | $r = 0.983$ |
| $SD(\text{md } 95) = 8.56$ | $Sy.x = 4.50$ |
- Number of samples measured: 49
The sample concentrations were between 58.6 and 216.3 mg/dl.

LDL-C Plus

References

- 1 Rifai N, Warnick GR, McNamara JR, Belcher JD, Grinstead GF, Frantz Jr ID. Measurement of Low-Density-Lipoprotein Cholesterol in Serum: a Status Report. Clin Chem 1992;38:150-160.
- 2 Cremer P, Seidel D. Lipoproteinanalytik: Methodische Empfehlungen. DG Klinische Chemie. Mitteilungen 1990;21:215-232.
- 3 National Cholesterol Education Program. Expert Panel on Detection, Evaluation, and Treatment of High Blood Cholesterol in Adults (Adult Treatment Panel II). NIH Publication No. 93-3095, 1995.
- 4 Naito HK, Strong JP, Scott MG, Roheim PS, Asztalos BF, Zilversmit DB, Srinivasan SR, Berenson GS, Wilson PWF, Scanu AM, Malikow MR. Atherogenesis: current topics on etiology and risk factors. Clin Chem 1995;41:132-133 No. 1.
- 5 Wieland H, Seidel D. Quantitative Lipoprotein Electrophoresis. In: Handbook of Electrophoresis, Vol III, ed. Lewis A., Boca Raton: CRC Press, 83-102, 1983.
- 6 Armstrong V, Seidel D. Evaluation of a Commercial Kit for the Determination of LDL-Cholesterol in Serum Based on Precipitation of LDL with Dextran Sulfate. Ärztl. Lab. 1985;31:325-330.
- 7 Friedewald WF, Levy RI, Frederickson DS. Estimation of LDL-Cholesterol Concentration without Use of the Preparative Ultracentrifuge. Clin Chem 1972;18:499-502.
- 8 Bachorik PS, Ross JW. National cholesterol education program recommendations for measurement of low-density lipoprotein cholesterol: executive summary. Clin Chem 1995;41:1414-1420.
- 9 Cohn JS, McNamara JR, Schaefer EJ. Lipoprotein Cholesterol Concentrations in the Plasma of Human Subjects as Measured in the Fed and Fasted States. Clin Chem 1988;34:2456-2459.
- 10 Pisani T, Gebiski CP, Leary ET, et al. Accurate Direct Determination of Low-Density Lipoprotein Cholesterol Using an Immuno-separation Reagent and Enzymatic Cholesterol Assay. Arch Pathol Lab Med 1995;119:1127.
- 11 Data on file at Roche.
- 12 LDL Cholesterol Method Certification Protocol for Manufacturers. National Reference System for Cholesterol. Cholesterol Reference Method Laboratory Network. October 1997.
- 13 Tietz, NW. Clinical Guide to Laboratory Tests. 3rd ed. WB Saunders Co, Philadelphia, Pa. 1995:404-407.
- 14 Cremer P, Nagel D, Mann H et al. Ten-year follow-up results from the Goettingen Risk, Incidence and Prevalence Study (GRIPS). I. Risk factors for myocardial infarction in a cohort of 5790 men. Atherosclerosis 1997;129:221-230.
- 15 Passing H, Bablok W. A New Biometrical Procedure for Testing the Equality of Measurements from Two Different Analytical Methods. J Clin Chem Clin Biochem 1983;21:709-720.
- 16 Bablok W et al. A General Regression Procedure for Method Transformation. J Clin Chem Clin Biochem 1988;26:783-790.

Instrument settings

US users

Refer to application sheet for additional operating information.

Roche/Hitachi 736 and 914 customers

Refer to application sheet for parameters.

Roche/Hitachi 904, 911, 912, 917 and MODULAR users

Read in the application parameters from the application diskette or barcode sheet, as appropriate.

Roche/Hitachi 704

Temperature: 37°C

PROGRAM 2 CHEMISTRY PARAMETERS	
TEST	[LDL]
ASSAY CODE	[1(1 POINT)] - [31] - [0]
SAMPLE VOLUME	[4]
R1 VOLUME	[300] - [50] - [NO]
R2 VOLUME	[100] - [20] - [NO]
WAVELENGTH	[700] - [600]
CALIB. METHOD	[LINEAR] - [0] - [0]
STD. (1) CONC.-POS.	[] - []
STD. (2) CONC.-POS.	[] - []
STD. (3) CONC.-POS.	[0] - [0]
STD. (4) CONC.-POS.	[0] - [0]
STD. (5) CONC.-POS.	[0] - [0]
STD. (6) CONC.-POS.	[0] - [0]
UNIT	[]
SD LIMIT	[0.1]
DUPLICATE LIMIT	[100]
SENSITIVITY LIMIT	[700]
ABS. LIMIT (INC/DEC)	[0] - [INCREASE]
PROZONE LIMIT	[32000] - [UPPER]
EXPECTED VALUE	[] - []
INSTRUMENT FACTOR	[1.00]


— Data entered by the operator

Roche/Hitachi 717

Temperature: 37°C

PROGRAM 2 CHEMISTRY PARAMETERS	
TEST	[LDL]
ASSAY CODE	[1(1 POINT)] - [50] - [0]
SAMPLE VOLUME	[4] - [3]
R1 VOLUME	[300] - [50] - [NO]
R2 VOLUME	[100] - [20] - [NO]
WAVELENGTH	[700] - [600]
CALIB. METHOD	[LINEAR] - [0] - [0]
STD. (1) CONC.-POS.	[] - []
STD. (2) CONC.-POS.	[] - []
STD. (3) CONC.-POS.	[0] - [0]
STD. (4) CONC.-POS.	[0] - [0]
STD. (5) CONC.-POS.	[0] - [0]
STD. (6) CONC.-POS.	[0] - [0]
SD LIMIT	[999]
DUPLICATE LIMIT	[100]
SENSITIVITY LIMIT	[700]
ABS. LIMIT (INC/DEC)	[0] - [INCREASE]
PROZONE LIMIT	[32000] - [UPPER]
EXPECTED VALUE	[] - []
PANIC VALUE	[] - []
INSTRUMENT FACTOR	[1.00]

— Data entered by the operator

LDL-C Plus

Roche/Hitachi 737

Temperature: 37°C

SYSTEM PARAMETERS CHAPTER 9.0. (CHEMISTRY)			
TEST NAME		LDL	
1.	ASSAY CODE	END - 20 - 00	
2.	SAMPLE VOLUME (µl)	4	
3.	R 1 VOLUME (µl)	300	
4.	R 2 VOLUME (µl)	100	
5.	WAVELENGTH 1	600 nm	
	WAVELENGTH 2	700 nm	
6.	COMPENSATE LIMIT	0.0	
7.	CALIBRATION		
	REQ. NO	CALIB. ID	CONC
	1) 01	D SALINE	0
	2) 02	CALIB.	Assigned value
	3)		
	4)		
	5)		
	6)		
	7)		
8.	EQUATION NO (1-5)	1	
9.	FACTOR (FIXED)	
10.	UNIT FACTOR	1.00	
11.	ABS. LIMIT (RATE)	0	
	INC/DEC	INC	

— Data entered by the operator

Roche/Hitachi 747

Temperature: 37°C

PROGRAM 4.2 CHEMISTRY PARAMETERS			
TEST (LDL)		
ASSAY CODE	1 - 50 - 00		
WAVELENGTH	700 (SUB) - 600 (MAIN)		
	SERUM	URINE	
SAMPLE VOLUME (µl)	4 - 3	[] - []	
EXPECTED VALUE	0.0 - 00.0	[] - []	
PANIC VALUE -	[] - []	
ABS. LIMIT (INC/DEC)	0 - INC	[] - []	
PROZONE LIMIT	32000 - UPPER	[] - []	
	R1	R2	
R1/R2 VOLUME (µl)	300	100	
R1/R2 DUMMY INTERVAL	0	0	
DILUTION VOLUME (µl)	0		
CALIB. METHOD	LINEAR		
POINTS	2		
STD 1 CONC RACK POS (NaCl)	0.0 - - 0		
STD 2 CONC RACK POS	Value - -		
STD 3 CONC RACK POS	0 - 0 - 0		
STD 4 CONC RACK POS	0 - 0 - 0		
STD 5 CONC RACK POS	0 - 0 - 0		
STD 6 CONC RACK POS	0 - 0 - 0		
SD LIMIT	999		
DUPLICATE LIMIT	100		
SENSITIVITY LIMIT	700		
STD 1 ABS. LEVEL	[-100]- [500].....		
INSTRUMENT FACTOR	1.00		

..... Data entered by the operator

Roche/Hitachi 902

No.	<Chemistry>	
1	Test Name	LDL
2	Assay Code (Mthd)	1 Point
3	Assay Code (2. Test)	0
4	Reaction Time	10
5	Assay Point 1	35
6	Assay Point 2	0
7	Assay Point 3	0
8	Assay Point 4	0
9	Wavelength (SUB)	700
10	Wavelength (MAIN)	600
11	Sample Volume	4.0
12	R1 Volume	300
13	R1 Position
14	R1 Bottle Size	Large
15	R2 Volume	0
16	R2 Position	0
17	R2 Bottle Size	Large
18	R3 Volume	100
19	R3 Position
20	R3 Bottle Size	Small
21	Calib. Type (Type)	Linear
22	Calib. Type (Wght)	2
23	Calib. Conc. 1	0
24	Calib. Pos. 1
25	Calib. Conc. 2
26	Calib. Pos. 2
27	Calib. Conc. 3	0
28	Calib. Pos. 3	0
29	Calib. Conc. 4	0
30	Calib. Pos. 4	0
31	Calib. Conc. 5	0
32	Calib. Pos. 5	0
33	Calib. Conc. 6	0
34	Calib. Pos. 6	0
35	S1 ABS.	0
36	K Factor	00000
37	K2 Factor	10000
38	K3 Factor	10000
39	K4 Factor	10000
40	K5 Factor	10000
41	A Factor	0
42	B Factor	0
43	C Factor	0
44	SD Limit	0.1
45	Duplicate Limit	100
46	Sens. Limit	700
47	S1ABS. Limit (L)	-100
48	S1ABS. Limit (H)	500
49	ABS. Limit	0
50	ABS. Limit (D/I)	Increase
51	Prozone Limit	32000
52	Proz Limit (Upp/Low)	Upper
53	Prozone (End Point)	00
54	Expect. Value (L)
55	Expect. Value (H)
56	Instr. Fact. (a)	1.0
57	Instr. Fact. (b)	0.0
58	Key Setting

.... Data entered by the operator

For detailed information, consult the operator manuals for Roche/Hitachi systems, the respective application sheets and the package inserts for the calibrators and control sera.

Tina-quant, Precinorm and Precipath are trademarks of a member of the Roche Group.

Intralipid is a trademark of KabiPharmacia, Inc.

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Bibliografia

- 1 Rifal N, Warnick GR, McNamara JR, Belcher JD, Grinstead GF, Frantz Jr ID. Measurement of Low-Density-Lipoprotein Cholesterol in Serum: a Status Report. Clin Chem 1992;38:150-160.
- 2 Cremer P, Seidel D. Lipoproteinanalytik: Methodische Empfehlungen. DG Klinische Chemie. Mitteilungen 1990;21:215-232.
- 3 National Cholesterol Education Program. Expert Panel on Detection, Evaluation, and Treatment of High Blood Cholesterol in Adults (Adult Treatment Panel II). NIH Publication No. 93-3095, 1995.
- 4 Naito HK, Strong JP, Scott MG, Roheim PS, Asztalos BF, Zilversmit DB, Srinivasan SR, Berenson GS, Wilson PWF, Scanu AM, Malikow MR. Atherogenesis: current topics on etiology and risk factors. Clin Chem 1995;41:132-133 No. 1.
- 5 Wieland H, Seidel D. Quantitative Lipoprotein Electrophoresis. In: Handbook of Electrophoresis, Vol III, ed. Lewis A., Boca Raton: CRC Press, 83-102, 1983.
- 6 Armstrong V, Seidel D. Evaluation of a Commercial Kit for the Determination of LDL-Cholesterol in Serum Based on Precipitation of LDL with Dextran Sulfate. Ärztl. Lab. 1985;31:325-330.
- 7 Friedewald WF, Levy RI, Frederickson DS. Estimation of LDL-Cholesterol Concentration without Use of the Preparative Ultracentrifuge. Clin Chem 1972;18:499-502.
- 8 Bachorik PS, Ross JW. National cholesterol education program recommendations for measurement of low-density lipoprotein cholesterol: executive summary. Clin Chem 1995;41:1414-1420.
- 9 Cohn JS, McNamara JR, Schaefer EJ. Lipoprotein Cholesterol Concentrations in the Plasma of Human Subjects as Measured in the Fed and Fasted States. Clin Chem 1988;34:2456-2459.
- 10 Pisani T, Gebiski CP, Leary ET, y cols. Accurate Direct Determination of Low-Density Lipoprotein Cholesterol Using an Immuno-separation Reagent and Enzymatic Cholesterol Assay. Arch Pathol Lab Med 1995;119:1127.
- 11 Documentação da Roche.
- 12 LDL Cholesterol Method Certification Protocol for Manufacturers. National Reference System for Cholesterol. Cholesterol Reference Method Laboratory Network. October 1997.
- 13 Tietz, NW. Clinical Guide to Laboratory Tests. 3ª ed. WB Saunders Co, Philadelphia, Pa. 1995:404-407.
- 14 Cremer P, Nagel D, Mann H et al. Ten-year follow-up results from the Goettingen Risk, Incidence and Prevalence Study (GRIPS). I. Risk factors for myocardial infarction in a cohort of 5790 men. Atherosclerosis 1997;129:221-230.
- 15 Passing H, Bablok W. A New Biometrical Procedure for Testing the Equality of Measurements from Two Different Analytical Methods. J Clin Chem Clin Biochem 1983;21:709-720.
- 16 Bablok W et al. A General Regression Procedure for Method Transformation. J Clin Chem Clin Biochem 1988;26:783-790.

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© 1999 Roche Diagnostics
Roche Diagnostics GmbH, D-68298 Mannheim, Alemanha