# $\alpha$ 1-Antitrypsin

Tina-quant®  $\alpha$ 1-Antitrypsin

● Indicates Roche/Hitachi analyzer(s) on which kit(s) can be used

Cat. no.	Bottle	Contents	704	717	736	747	902	904	911	914	917	MODULAR	
					737							P	D
1557599	1	Buffer, 2 x 15 ml	●	●			●	●	●	●	●	●	
	2	$\alpha$ 1-Antitrypsin antibodies, 2 x 6 ml											

Some analyzers and kits shown may not be available in all countries. For additional system applications, contact your local Roche representative.

## Intended use

Immunoturbidimetric assay for the quantitative in vitro determination of  $\alpha$ 1-antitrypsin in human serum and plasma on automated clinical chemistry analyzers.

## Summary<sup>1-6</sup>

$\alpha$ 1-Antitrypsin is a glycoprotein consisting of a polypeptide chain to which 3 oligosaccharide chains are bonded (54000 daltons). It is synthesized in hepatocytes.  $\alpha$ 1-Antitrypsin is quantitatively the most important proteinase inhibitor (Pi) in serum and plasma. It specifically inactivates serine proteases (e.g. trypsin, chymotrypsin, collagenase, leukocyte elastase, plasmin and thrombin), with which it reversibly forms an enzyme-indicator complex.  $\alpha$ 1-Antitrypsin constitutes the major fraction of the electrophoretic  $\alpha$ 1-globulin fraction. Because it is a small molecule,  $\alpha$ 1-antitrypsin diffuses rapidly into other body fluids including bronchial secretions.

$\alpha$ 1-Antitrypsin is an important, positive acute phase reactant found in elevated concentrations in inflammatory processes (e.g. infectious and rheumatoid diseases), tissue necrosis, malignancy and traumas. Inflammation of the liver parenchymal cells is often accompanied by elevated  $\alpha$ 1-antitrypsin levels, although levels of other acute phase reactants are not affected.


Acute hereditary  $\alpha$ 1-antitrypsin deficiency is suspected in cases of neonatal hepatitis accompanied by progressive liver cirrhosis in early childhood. It is also suspected when severe pulmonary emphysema occurs in adults due to the prevalence of leukocyte elastase, which can lead to unrestrained proteolytic degradation of the pulmonary parenchymal cells.

Various methods are available for assaying  $\alpha$ 1-antitrypsin, such as kinetic nephelometry or radial immunodiffusion. Roche's automated  $\alpha$ 1-antitrypsin assay is based upon the immunological agglutination test principle.

## Test principle<sup>1</sup>

Immunoturbidimetric assay

- Sample and addition of R1 (buffer)
- Addition of R2 ( $\alpha$ 1-antitrypsin antibody) and start of reaction:

Antibodies to  $\alpha$ 1-antitrypsin react with antigen in the sample to form an antigen/antibody complex which can be determined turbidimetrically after agglutination.

## Working solution concentration

### R1 Buffer

Phosphate buffer: 12.7 mmol/l, pH 7.2; NaCl: 0.13 mol/l; PEG: 50 g/l; preservative

### R2 $\alpha$ 1-Antitrypsin antibody solution

Anti-human  $\alpha$ 1-antitrypsin antibody (rabbit) > 2 g/l; NaCl: 0.12 mol/l; preservative

## Precautions and warnings

For in vitro diagnostic use.

Exercise the normal precautions required for handling all laboratory reagents.

## Reagent handling

- R1: Ready for use
- R2: Ready for use

## Storage and stability

Unopened kit components: Up to the expiration date at 2–8°C

R1: 90 days opened and refrigerated on the analyzer

R2: 90 days opened and refrigerated on the analyzer

## Specimen collection and preparation

Collect serum using standard sampling tubes.

Heparinized, EDTA or citrated plasma

Stability<sup>7</sup>: 7 days at 20–25°C

3 months at 4–8°C

3 months at -20°C

Thaw only once.

Roche/Hitachi analyzers with automatic sample dilution

Serum or plasma samples are diluted **1 + 15** on the analyzer with 0.9% NaCl.

Roche/Hitachi analyzers without automatic sample dilution

Dilute serum or plasma samples manually **1 + 15** with 0.9% NaCl.

$\alpha$ 1-Antitrypsin, ceruloplasmin, haptoglobin and prealbumin can be determined selectively in the diluted sample.

Centrifuge samples containing precipitate before performing the assay.

## Testing procedure

### Materials provided

- Working solutions as described above

### Additional materials required

- Calibrators and controls as indicated below
- 0.9% NaCl

## Assay

Refer to the appropriate operator's manual and/or Instrument Settings section of this insert for analyzer-specific assay instructions. The performance of applications not validated by Roche is not warranted and must be defined by the user.

## Calibration

Standardization: The  $\alpha$ 1-antitrypsin method was standardized against the reference preparation CRM 470 corresponding to RPPHS (Reference Preparation for Proteins in Human Serum).

Roche/Hitachi 704/717/902/914

S1: 0.9% NaCl

S2–6: Diluted C.f.a.s (Calibrator for automated systems) Proteins

The calibration is performed using C.f.a.s Proteins via a dilution series (6-point calibration).

Dilute the C.f.a.s. Proteins calibrator in the following manner to prepare S2–6:

Std. No.	Dilution Sample + 0.9% NaCl	Conversion factor
2	1 + 127	0.125
3	1 + 63	0.250
4	1 + 31	0.500
5	1 + 15	1.000
6	1 + 6	2.286

Multiply the lot-specific C.f.a.s Proteins calibrator value by the factors above to determine the standard concentrations for the 6-point calibration curve.

# $\alpha$ 1-Antitrypsin

Roche/Hitachi 904/911/912/917/MODULAR

S1: 0.9% NaCl

S2: C.f.a.s. Proteins

Multiply the lot-specific C.f.a.s. Proteins calibrator value by the factors below to determine the standard concentrations for the six-point calibration curve.

S2: 0.125            S5: 1.000

S3: 0.250            S6: 2.286

S4: 0.500

Calibration frequency

Full calibration is recommended:

- after lot change
- as required following quality control procedures

Calibration verification: Not necessary.

## Quality control

For quality control, use Precinorm Protein and Precipath Protein or other suitable control material.

Roche/Hitachi analyzers with automatic sample dilution:

The control material is diluted **1 + 15** on the analyzer with 0.9% NaCl in the same way as the samples.

Roche/Hitachi analyzers without automatic sample dilution:

Dilute the control material manually **1 + 15** with 0.9% NaCl in the same way as the samples.

The control intervals and limits must be adapted to the individual laboratory and country-specific requirements. Values obtained should fall within established limits. Each laboratory should establish corrective measures to be taken if values fall outside the limits.

## Calculation

Roche/Hitachi systems automatically calculate the  $\alpha$ 1-antitrypsin concentration of each sample.

Conversion factor: mg/dl  $\times$  0.01 = g/l

g/l  $\times$  100 = mg/dl

## Limitations – interference<sup>8</sup>

Criterion: Recovery within  $\pm$  10% of initial values.

Icterus: No significant interference up to an I index of 60 (approximate conjugated and unconjugated bilirubin concentration: 60 mg/dl).

Hemolysis: No significant interference up to an H index of 1000 (approximate hemoglobin concentration: 1000 mg/dl).

Lipemia (Intralipid): No significant interference up to an L index of 890 (approximate triglycerides concentration: 1780 mg/dl). There is poor correlation between turbidity and triglycerides concentration.

Rheumatoid factors < 100 IU/ml do not interfere.

No high dose hook effect is seen up to an  $\alpha$ 1-antitrypsin concentration of 1200 mg/dl.

Elevated estrogen levels (oral contraceptives; third trimester of pregnancy) give rise to false positive results. CRP and haptoglobin should therefore also be determined.<sup>9</sup>

## Measuring/reportable range

Roche/Hitachi 704

Measuring range: 30–600 mg/dl (0.30–6.00 g/l)\*

Roche/Hitachi 902

Measuring range: 30–420 mg/dl (0.30–4.20 g/l)\*

Roche/Hitachi 717/904/911/912/914/917/MODULAR

Measuring range: 30–600 mg/dl (0.30–6.00 g/l)\*

Extended measuring range with rerun: 30–1200 mg/dl (0.30–12.0 g/l)\*\*

Determine samples with higher concentrations via the rerun function. On instruments without rerun function, manually dilute samples with higher concentrations using 0.9% NaCl (e.g. 1 + 1). Multiply the result by the appropriate dilution factor (e.g. 2).

\* Maximum reportable range is dependent on the highest standard concentration.

\*\* Maximum extended reportable range is an approximate value dependent on the calculated value of the highest standard.

## Expected values\*\*\*

90–200 mg/dl (0.9–2.0 g/l)

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference ranges. For diagnostic purposes,  $\alpha$ 1-antitrypsin results should always be assessed in conjunction with the patient's medical history, clinical examinations and other findings.

\*\*\* Reference range in accordance with CRM 470 protein standardization<sup>10</sup>

## Specific performance data

The data determined using a Roche/Hitachi system are given below. The results obtained in individual laboratories may differ.

## Imprecision<sup>8</sup>

Reproducibility was determined using human samples and controls in an internal protocol: n = 21. The following results were obtained.

Sample	Within run			Between day		
	Mean mg/dl	SD mg/dl	%CV	Mean mg/dl	SD mg/dl	%CV
Human serum	160.0	1.8	1.1	154.6	4.7	3.0
Precinorm Protein	154.9	1.1	0.7	156.3	5.7	3.6
Precipath Protein	307.8	4.7	1.5	314.7	12.2	3.9

## Analytical sensitivity (lower detection limit)<sup>8</sup>

Detection limit: 10 mg/dl (0.10 g/l)

The lower detection limit represents the lowest measurable concentration of  $\alpha$ 1-antitrypsin that can be distinguished from zero. It is calculated as three standard deviations of 21 replicates of the lowest standard.

## Analytical specificity<sup>8</sup>

The  $\alpha$ 1-antitrypsin antibody is specific for the  $\alpha$ 1-antitrypsin protein.

## Method comparison<sup>8</sup>

A comparison of the  $\alpha$ 1-antitrypsin determination using the Roche Tina-quant  $\alpha$ 1-antitrypsin assay (y) with a nephelometric  $\alpha$ 1-antitrypsin assay (x) in human sera gave the following correlation (mg/dl):

Passing/Bablok<sup>11,12</sup>            Linear regression

y = 9.9 + 0.993 x            y = 12.2 + 0.971 x

r = 0.967                        r = 0.967

SD (md 95) = 20.2            Sy.x = 11.4

Number of samples measured: 123

The sample concentrations were between 91.5 and 471.7 mg/dl.

## References

- 1 Lievens M, Bienvenu J, Buitrago JMG et al. Evaluation of four new Tina-quant  $\alpha$ 1-antitrypsin assays for determination of  $\alpha$ 1-acid glycoprotein,  $\alpha$ 1-antitrypsin, haptoglobin and prealbumin. Clin Lab 1996;42: 515–520.
- 2 Greiling H, Gressner AM, eds. Lehrbuch der klinischen Chemie und Pathobiochemie, 3rd ed. Stuttgart/New York: Schattauer, 1995:231–232.
- 3 Schneider M, Pott G, Gerlach U. Alpha-1-Antitrypsin-Mangel. Klin Wschr 1986;64:197–205.
- 4 Tietz NW, ed. Clinical Guide to Laboratory Tests, 3rd ed. Philadelphia, Pa: WB Saunders, 1995:66–67.
- 5 Johnson AM, Alper CA. Deficiency of  $\alpha$ 1-antitrypsin in childhood liver disease. Pediatrics 1970;46:921–925.
- 6 Silverman EK, Pierce JA, Province MA et al. Variability of pulmonary function in alpha-1-antitrypsin deficiency: Clinical correlation. Ann Intern Med 1989;111:982–991.
- 7 Guder WG, Narayanan S, Wissner H, Zawta B. List of Analytes; Pre-analytical Variables. In: Samples: From the Patient to the Laboratory. Darmstadt: GIT-Verlag, 1996.
- 8 Data on file at Roche.
- 9 Laurell CB, Kullander S, Thorell J. Effect of administration of a combined estrogen-progestin contraceptive on the level of individual plasma proteins. Scand J Lab Invest 1968;21:337–343.
- 10 Consensus values of the "Deutsche Gesellschaft für Laboratoriumsmedizin", the "Deutsche Gesellschaft für Klinische Chemie" and the "Verband der Diagnostica-Industrie e.V." (VDGH). DG Klinische Chemie Mitteilungen 1995;26:119–122.
- 11 Passing H, Bablok W. A New Biometrical Procedure for Testing the Equality of Measurements from Two Different Analytical Methods. J Clin Chem Clin Biochem 1983;21:709–720.
- 12 Bablok W et al. A General Regression Procedure for Method Transformation. J Clin Chem Clin Biochem 1988;26:783–790.

# $\alpha$ 1-Antitrypsin

English - 1999-09 - 2133717001 03 02

## Instrument settings

### US users

Refer to application sheet for additional operating information.

### Roche/Hitachi 914 customers

Refer to application sheet for parameters.

### Roche/Hitachi 904, 911, 912, 917, and MODULAR users

Read in the application parameters from the application diskette or barcode sheet, as appropriate.

### Roche/Hitachi 704

Temperature: 37°C

PROGRAM 2 CHEMISTRY PARAMETERS	
TEST	[ATRYP]
ASSAY CODE	[2(2 point)]-[15]-[32]
SAMPLE VOLUME	[14]
R1 VOLUME	[340]-[20]-[NO]
R2 VOLUME	[113]-[20]-[NO]
WAVELENGTH	[700]-[340]
CALIB. METHOD	[NON-LINEAR]-[1]-[6]
STD. (1) CONC.-POS.	[ ]-[ ]
STD. (2) CONC.-POS.	[ ]-[ ]
STD. (3) CONC.-POS.	[ ]-[ ]
STD. (4) CONC.-POS.	[ ]-[ ]
STD. (5) CONC.-POS.	[ ]-[ ]
STD. (6) CONC.-POS.	[ ]-[ ]
UNIT	[ ]
SD LIMIT	[100]
DUPLICATE LIMIT	[300]
SENSITIVITY LIMIT	[0]
ABS. LIMIT (INC/DEC)	[0]-[INCREASE]
PROZONE LIMIT	[32 000]-[UPPER]
EXPECTED VALUE	[ ]-[ ]
INSTRUMENT FACTOR	[1.00]

\_\_\_ Data entered by the operator

### Roche/Hitachi 717

Temperature: 37°C

PROGRAM 2 CHEMISTRY PARAMETERS	
TEST	[ATRYP]
ASSAY CODE	[2(2 point)]-[24]-[50]
SAMPLE VOLUME	[10]-[5]
R1 VOLUME	[250]-[20]-[NO]
R2 VOLUME	[83]-[20]-[NO]
WAVELENGTH	[700]-[340]
CALIB. METHOD	[NON-LINEAR]-[1]-[6]
STD. (1) CONC.-POS.	[ ]-[ ]
STD. (2) CONC.-POS.	[ ]-[ ]
STD. (3) CONC.-POS.	[ ]-[ ]
STD. (4) CONC.-POS.	[ ]-[ ]
STD. (5) CONC.-POS.	[ ]-[ ]
STD. (6) CONC.-POS.	[ ]-[ ]
SD LIMIT	[100]
DUPLICATE LIMIT	[300]
SENSITIVITY LIMIT	[0]
ABS. LIMIT (INC/DEC)	[0]-[INCREASE]
PROZONE LIMIT	[32 000]-[UPPER]
EXPECTED VALUE	[ ]-[ ]
PANIC VALUE	[ ]-[ ]
INSTRUMENT FACTOR	[1.00]

\_\_\_ Data entered by the operator

## Roche/Hitachi 902

No.	<Chemistry>	
1	Test Name	ATRYP
2	Assay Code (Mthd)	2 Point End
3	Assay Code (2. Test)	0
4	Reaction Time	10
5	Assay Point 1	17
6	Assay Point 2	35
7	Assay Point 3	0
8	Assay Point 4	0
9	Wavelength (SUB)	700
10	Wavelength (MAIN)	340
11	Sample Volume	10.0
12	R1 Volume	250
13	R1 Pos.	.....
14	R1 Bottle Size	Small
15	R2 Volume	0
16	R2 Pos.	0
17	R2 Bottle Size	Small
18	R3 Volume	83
19	R3 Pos.	.....
20	R3 Bottle Size	Small
21	Calib. Type (Type)	Logit-log(4P)
22	Calib. Type (Wght)	0
23	Calib. Conc. 1	0.0
24	Calib. Pos. 1	.....
25	Calib. Conc. 2	.....
26	Calib. Pos. 2	.....
27	Calib. Conc. 3	.....
28	Calib. Pos. 3	.....
29	Calib. Conc. 4	.....
30	Calib. Pos. 4	.....
31	Calib. Conc. 5	.....
32	Calib. Pos. 5	.....
33	Calib. Conc. 6	.....
34	Calib. Pos. 6	.....
35	S1 ABS.	0
36	K Factor	10000
37	K2 Factor	10000
38	K3 Factor	10000
39	K4 Factor	10000
40	K5 Factor	10000
41	A Factor	0
42	B Factor	0
43	C Factor	0
44	SD Limit	100
45	Duplicate Limit	300
46	Sens. Limit	0
47	S1ABS. Limit (L)	-32000
48	S1ABS. Limit (H)	32000
49	ABS. Limit	32000
50	ABS. Limit (D/L)	Increase
51	Prozone Limit	32000
52	Proz Limit (Upp/Low)	Upper
53	Prozone (End Point)	35
54	Expect. Value (L)	.....
55	Expect. Value (H)	.....
56	Instr. Fact. (a)	1.0
57	Instr. Fact. (b)	0.0
58	Key Setting	.....

..... Data entered by the operator

# $\alpha$ 1-Antitrypsin

For detailed information, consult the operator manuals for Roche/Hitachi systems, the respective application sheets and the package inserts for the calibrator and control sera. Tina-quant, Precinorm and Precipath are trademarks of a member of the Roche Group. Intralipid is a trademark of KabiPharmacia, Inc.

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© 1999 Roche Diagnostics
Roche Diagnostics GmbH, D-68298 Mannheim