Creatine Kinase

Liquid Reagent

Creatine Quinase

Reagente líquido

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COBAS INTEGRA Creatine Kinase	200 testes	Ref. 04524977 190 System-ID 07 5923 6
Calibrator f.a.s.	12 × 3 mL	Ref. 10759350 190
Calibrator f.a.s. (para EUA)	12 × 3 mL	Ref. 10759350 360 System-ID 07 3718 6
Precinorm U	20 × 5 mL	Ref. 10171743 122 System-ID 07 7997 0
Precipath U	20 × 5 mL	Ref. 10171778 122 System-ID 07 7998 9
Precinorm U plus	10 × 3 mL	Ref. 12149435 122 System-ID 07 7999 7
Precipath U plus	10 × 3 mL	Ref. 12149443 122 System-ID 07 8000 6

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Bibliografia

- Lott JA, Stang JM. Serum enzymes and isoenzymes in the diagnosis and differential diagnosis of myocardial ischemia and necrosis. Clin Chem. 1980;26:1241-1250.
- Moss DW, Henderson AR, Kachmar JF. Enzymes. In: Tietz NW, ed. Fundamentals of Clinical Chemistry. 3rd ed. Philadelphia: WB Saunders 1987:346-421.
- Hörder M, Elser RC, Gerhardt W, Mathieu M, Sampson EJ. IFCC methods for the measurement of catalytic concentration of enzymes. Part 7. IFCC method for creatine kinase (ATP: creatine N-phosphotransferase, EC 2.7.3.2). J Int Fed Clin Chem 1989;1:130-139.
- Mathieu M, Breaudiere JP, Galteau MM, Giudollet J, Lalegerie P, Bailly M, et al. Recommendations for measuring the catalytic concentration of creatine kinase in human serum at 30°C. Ann Biol Clin 1982;40:87-164.
- Hörder M, Magid E, Pitkänen E, Härkönen M, Strömme JH, Theodorsen L, et al. Recommended method for the determination of creatine kinase in blood modified by the inclusion of EDTA. Scand J Clin Lab Invest 1979;39:1-5.
- Bergmeyer HU, Breuer H, Büttner H, Delbrück A, Laue D, Pilz W, et al. Empfehlungen der Deutschen Gesellschaft für Klinische Chemie. Standard-Methode zur Bestimmung der Aktivität der Creatin-Kinase. J Clin Chem Clin Biochem 1977;15:249-254.
- Guder WG, Narayanan S, Wissner H, Zawta B. List of Analytes; Pre-analytical Variables. Brochure in: Samples: From the Patient to the Laboratory. Darmstadt: GIT Verlag 1996.
- Glick MR, Ryder KW, Jackson SA. Graphical Comparisons of Interferences in Clinical Chemistry Instrumentation. Clin Chem 1986;32:470-474.
- Breuer J, Report on the Symposium "Drug Effects in Clinical Chemistry Methods". Eur J Clin Chem Clin Biochem 1996;34:385-386.

10. Klein G, Berger A, Bertholf R et al. Abstract: Multicenter Evaluation of Liquid Reagents for CK, CK-MB and LDH with Determination of Reference Intervals on Hitachi Systems. Clin Chem 2001;47:Suppl. A30.
11. Stein W. Strategie der klinischen-chemischen Diagnostik des frischen Myokardinfarktes. Med Welt 1985;36:572-577.
12. Thomas L, Müller M, Schumann G, Weidemann G et al. Consensus of DGKL and VDGH for interim reference intervals on enzymes in serum. J Lab Med 2005;29:301-308.
13. Tietz NW. Clinical Guide to Laboratory Tests, 3rd edition. Philadelphia, PA: WB Saunders Co, 1995:180-181.
14. Myocardial Infarction Redefined - A Consensus Document of the Joint European Society of Cardiology/American College of Cardiology Committee for the Redefinition of Myocardial Infarction. Eur Heart J 2000;21:1502-1513.
15. Black HR, Quallich H, Gareleck CB. Racial differences in serum creatine kinase levels. Am J Med 1986;81:479-487.
16. Bablok W et al. A General Regression Procedure for Method Transformation. J Clin Chem Clin Biochem 1988;26:783-790.

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©2006 Roche Diagnostics.
Roche Diagnostics GmbH, D-68298 Mannheim